

PIÈCES DÉTACHÉES A-DEC 2015 Q1 Authorized Dealer Catalog





A-dec Catalogue des pièces détachées 2015

Seules des pièces d'origine A-dec peuvent maintenir
votre équipement A-dec comme neuf.



Table des matières

INFORMATION GÉNÉRALE

Garantie pour les produits A-dec

PÉRIPHÉRIQUES

Pièces pour moteurs électriques.....	1
Source lumière, ampoules	2
Réglage air pièce à main.....	3

EAU PROPRE AUTONOME

A-dec 500 et 300	4
Cascade et avant.....	7

COMPOSANTS

Diaphragmes et joints.....	8
Joints toriques.....	9
Caches et couvercles A-dec 500	11
Caches et couvercles A-dec 300	12

EMBOUTS

Embout et manchons.....	14
Attaches rapides.....	17

BOÎTES AU SOL

Profilées.....	18
Inoxidables	19
Bouchons et tuyau ombilical.....	20
Séparateurs d'eau et filtres	21
Filtres/régulateurs.....	22
Filtres, pré-régulateurs, manomètres, valve d'arrêt	23

OUTILS

Clés hexagonales, test anti-rétraction.....	24
Peinture de retouche et silicone	25

TUBULURES ET TUYAUX

Codification des tubulures 1/8°	26
Codification des tubulures suite	27
Cordons en silicone.....	28
Cordons en vinyl	29
Cordons en vinyl, raccords et ampoules.....	30
A-dec 500 et 300 Aspiration grand volume.....	31
A-dec 500 et 300 Pompe à salive.....	32
Cascade et Performer Aspiration grosse canule.....	33

A-dec 500, 300 et cascade valves d'aspiration 15 mm.....	34
Tuyaux d'aspiration, pédale au sol, boîte au sol et crachoir	35

VALVES ET CLAPETS

Momentanées.....	36
Manuel ON/OFF, réglage spray	37
Robinet de réglage	38
Relais	39
Coupure d'eau	40
Valves de sélection.....	41
Étrangleurs, leviers en plastique	42

SERINGUES

Canules et embouts seringues	43
Autoclavables traditionnelles	44
Autoclavables continentales.....	45
Boutons et kits de réparation	46
Kits de réparation.....	47
Composants et tuyaux	48

PÉDALES PNEUMATIQUES

Pédales montées.....	49
Composants des pédales et kits de reconstruction	50

CARTOUCHES, CLAPETS ET SUPPORTS

A-dec 500 & 300 cartouches, clapets, diaphragmes	53
Bloc Century plus	54
Bloc Century II.....	56
Bloc Performer.....	58
Fouets, poignées et supports.....	59
Supports pièces à main et collecteur	60
Valves de pré-sélection	61

COMPOSANTS ÉLECTRIQUES

A-dec 500 et 300	62
A-dec 500 et 300 câbles RJ, électrovanne et module relais.....	63
Platines électriques, 25 W, 80 W et 300 W	64
Câbles d'alimentation pour 541, 542 et 545	65
Claviers	66
Interrupteurs électriques	67

COMPOSANTS ASPIRATION

Pompes à salive et grosse canule, adaptateur	68
--	----

Écouvillons, entonnoir, canules chir., kits de conversion.....	69
A-dec 500.....	70
A-dec 300.....	71
Collecteurs	72

CRACHOIRS

A-dec 500, 400 et 300.....	73
Cascade.....	74
Cascade.....	75
Performer	76
Composants	77

ÉCLAIRAGES

Composants de la LED 570	79
Composants de l'éclairage 300.....	80
A-dec 500, cascade, 6300 et Performer	81
Composants de l'éclairage 6300, transfo.	82
Câbles, transfo., interrupteurs et barres.....	83
Composants pour Performer.....	84

FAUTEUILS ET TABOURETS

Protège-pieds	85
Selleries vis et molettes	86
Remplacement des selleries.....	87
Vérins.....	90
Électrovannes et huile	91
Condensateurs, alimentations et platines.....	92
Capteurs de position	93
Potentiomètre et butée mécanique.....	94
Freins de rotation, câbles électriques et cordons d'alimentation.....	95
Pédales électriques.....	96
Barres de montage et accoudoirs	97
Molettes et roulettes avant 2015	99
Vérins et dossiers avant 2015.....	100

ACCESSOIRES

Supports plateaux	101
Plateaux et tapis.....	103
Protections units.....	104

INDEX

Index des références.....	105
Index alphabétique.....	111

Garantie pour les produits A-dec

Portée

A-dec garantit tous les produits répertoriés ci-dessous contre tout défaut de matériel ou de fabrication dans des conditions normales d'utilisation, achetés auprès d'A-dec ou d'un revendeur A-dec agréé. La seule obligation d'A-dec dans le cadre de cette garantie consiste à fournir des pièces pour la réparation des composants défectueux, ou leur remplacement par des composants équivalents en termes de performances. La présente garantie remplace toute autre garantie ou obligation, expresse ou implicite. A-dec exclut expressément toute garantie implicite de qualité marchande ou d'adéquation à un usage particulier. A-dec ne pourra en aucun cas être tenue pour responsable des dommages indirects, particuliers, accessoires, exemplaires ou consécutifs d'aucune sorte.

La réparation des composants ou leur remplacement dépendent des disponibilités en termes de couleurs et de styles au moment du recours en garantie. Si un recours en garantie est déposé après l'arrêt d'un produit ou d'une couleur, A-dec se réserve le droit d'honorer la garantie de l'une des manières suivantes :

1. Remplacer le composant concerné par un autre du même style et/ou de la même couleur
2. Remplacer le composant concerné et tout autre composant du produit par un autre de couleur et/ou de style similaires

L'acheteur n'aura aucun autre recours. Les garanties définies dans les présentes prévoient des droits spécifiques. D'autres droits peuvent exister, selon la juridiction d'achat et d'utilisation du produit, dans la mesure où ces droits s'appliquent au présent accord. Les exclusions de garantie sont fournies dans les sections suivantes. Les avis relatifs à la défaillance d'un produit doivent être adressés à A-dec durant la période de garantie.

Période

La période de garantie, qui court à partir de la date de livraison, est la suivante :

PRODUITS	GARANTIE
Cylindre hydraulique du fauteuil (inclinaison et levage)	10 ans
Fauteuils, tabourets de praticien et d'assistant, éclairages, unités, supports d'écran, mobilier dentaire et composants connexes	5 ans
Module de commande et moteurs électriques NLX d'A-dec (couverture selon la clause V. de la section Exclusions)	1 an
Cordons, modules de commande et moteur électrique EA-53 d'A-dec, Produits A-dec W&H Synea, moteurs pneumatiques et Assistina (couverture selon la clause V. de la section Exclusions)	2 ans
Produits A-dec W&H Alegria (couverture selon la clause V. de la section Exclusions)	1 an

Les composants remplacés sous garantie ou achetés conservent la garantie du produit d'origine ou sont garantis pendant une période d'un an après leur livraison (la période la plus longue l'emportant).

Exclusions

La garantie ne couvre pas les éléments suivants :

- I. Les dommages résultant d'une installation ou d'un entretien inadéquats, d'un accident, d'une mauvaise utilisation, d'un abus, d'une négligence, de modifications, du transport ou d'une catastrophe naturelle.
- II. Les dommages résultant de l'utilisation de produits chimiques et du processus de nettoyage, de désinfection ou de stérilisation.
- III. Les composants d'utilisation ordinaire tels que (mais sans s'y limiter) : écrans protecteurs de lampe, ampoules, filtres, joints, cordons et cartouches d'eau. De même, la garantie ne couvre pas l'usure normale.
- IV. Les produits consommables utilisés après la date d'expiration indiquée sur l'emballage.
- V. Tous les autres produits non fabriqués par A-dec, OEM inclus, sont couverts par la garantie du fabricant d'ori-

gine et non par celle d'A-dec. Exemples (mais sans s'y limiter) : stérilisateurs, matériel d'entretien, appareils photo, lampes à polymériser, instruments à ultrasons, modules de commande, moteurs électriques, fixations, instruments et turbines. Pour plus d'informations sur la garantie, contactez le service clientèle d'A-dec.

Le non-respect des instructions fournies dans le mode d'emploi A-dec (instructions d'utilisation et d'entretien) entraîne l'annulation de la garantie.

ATTENTION

La loi fédérale des États-Unis stipule que ce dispositif ne peut être vendu que par ou sur l'ordonnance d'un dentiste, d'un médecin ou d'un praticien diplômé en vertu de la loi de l'État dans lequel il exerce, et utilise ou prescrit l'utilisation de l'appareil.

Pièces pour moteurs électriques

EA-53 REPLACEMENT PARTS	PART NUMBER	PRICE
Motor Only, EA-53	WH3017900	
CM-2.2 module for up to 2 motors (handpiece tubing not included) Use with: EA-53	54.0480.01	
A-dec 82" VE-10 tubing Use with: EA-53	WH3018001	
Joint nez moteur Use with: EA-53	0.01862301	

EA-50, EA-51LT, EA-52LED REPLACEMENT PARTS	PART NUMBER	PRICE
Control Module, CM-2.0E Use with: EA-50LT, EA-51LT, EA-52LED	54.0451.01	
Tubing, VE-9, 82", Surf 6 Use with: EA-40LT, EA-50, EA-51LT, EA-52LED	WH04481100	
Tubing, VE-9, 82", Gray 3 Use with: EA-50, EA-51LT, EA-52LED	WH05343700	
Halogen Bulb Use with: EA-50, EA-51	WH06307200	
Upgrade Kit to LED Use with: EA-51LT, EA-50LT	WH06307200	
LED Bulb Use with: EA-52LED, upgraded, EA-50LT, EA-51LT	WH06189400	
Cordset O-ring Use with: EA-40LT, EA-50LT, EA-51LT, EA-52LED	WH01218700	
Platine LED pour EA52 Use with: EA-40LT, EA-50LT, EA-51LT, EA-52LED	WH06286500	



Source lumière, ampoules

Quad Voltage Intraoral Light Source




IOLS kit for A-dec 300


- Fiber-optic intraoral light source
- Requires fiber-optic handpiece tubing, 24vac
- All IOLS come equipped to light (2) handpieces
- A-dec 300 kit equipped to light (4) handpieces
- To add additional handpieces order the QVIOLS Switch Kit

PRODUCT DESCRIPTION	PART NUMBER	PRICE
Quad Volt IOLS kit for A-dec 300 (equipped to light (4) handpieces)	90.1168.00	

End Cap Parts




A-dec | W&H replacement bulbs




For 896, 898, 898-LE,
TA-96L, TA-98L



EA-40LT, EA-50LT,
EA-51LT, and A-25LT
electric motor



For LT25 lowspeed
motor and 6-pin to
5-hole adapter




Spacer for 6-pin
to 5-hole adapter

BULB BY HANDPIECE MODEL	PART NUMBER	PRICE
for EA-40LT, EA-50LT, EA-51LT, AM-25LT and A-25LT electric motor	0.04300400	

Réglage air pièce à main

A-dec Handpiece Exhaust Muffler


98.0660.00



- Attaches to any 4-hole vinyl handpiece tubing
- Particularly effective in quieting handpiece exhaust inside the control head


Handpiece Pressure Gauge (single)

50.0271.00



- Fit all handpiece couplers, except Borden 3-hole


A-dec 500 et 300



REF	DESCRIPTION	WHITE	SURF	PRICE
1	Water Bottle Assembly	14.0468.01	14.0468.00	
2	O-ring, AS568-006, 0.114 ID x 0.070 W	030.006.02		
3	Connector, water	14.0457.00		
4	Receptacle	14.0453.00		
5	Retaining ring, external tooth, .187" ID	010.002.00		
6	O-ring, N, 1.612 ID x .103 W (10)	031.130.01		
7	Tube, supply, 8" long	14.0467.00		
8	Clamp, water bottle	14.0458.01	14.0458.00	
9	Screw, button head socket, 10-32 x 1/2	005.110.00		
10	Water bottle with base	14.0466.01	14.0466.00	

Note: Use care when servicing this water bottle, be sure to align the clamp with the detent on the bottle to successfully reassemble the parts without causing damage.

A-dec 500 et 300



REF	DESCRIPTION	WHITE	SURF	PRICE
1-6	Kit, Remote Mount, Water Bottle includes 20' tubing	14.0471.00	-	
2-4, 6	Kit, Chair Mount, Water Bottle includes 9.5' tubing	14.0484.00	14.0464.00	
4, 6	Kit, Support Center / Slimline, Water Bottle includes 3' tubing	14.0464.02	-	
1	Casting bracket bracket mounts to flat surfaces	77.0420.01		
2	Housing, bottle	14.0455.04	14.0455.02	
3	Screw, button head socket, 10-32 x 3/8 stainless steel	001.253.00		
4	Bottle assembly, water	14.0468.01	14.0468.00	
5	Screw	003.071.00		
6	Receptacle assembly	90.1101.00		
7	Socket Head Screw 10-32 x 1/2	002.135.00		
8	O-ring, E, 0.921 ED x 0.139 W	034.213.00		
	Jacket tubing 5/16 x 1/8 (per foot)	024.200.00	-	

A-dec 500 et 300

A-dec ICX® Waterline 2-Liter Treatment Tablet

90.1065.00 Case, 2-liter tablets



- Case contains (36) boxes with (50) tablets each

Water Bottle Tubing Grommet

018.094.00



- For routing water bottle tubing through dental furniture
- 3/4" (19 mm)
- White

2-liter Water Bottle Cap

14.0491.01



Hole Plug for Unit

PART NUMBER

PRICE

Plug, no water bottle, Surf 4

77.0263.00

Plug, no water bottle, White 2

77.0263.01



- Plugs hole in front mount rigid arm when self-contained water bottle is not mounted

Water Bottle Brush

049.070.00



Cascade et avant

Self-Contained Water Systems with 0.7L Bottle




- Kits Include:
Regulator assembly
(2) 0.7L water bottles with tamper-proof caps
Bottle receptacle assembly
Mounting hardware
- Can be used with ICX or an approved maintenance protocol to maintain dental water lines
- Provides water for handpieces and syringe

14.0424.00 Radius 2122/2132 L/R Systems

Includes short stud
with 14.0492.00 style receptacle assembly

Replacement Parts for Self-Contained 0.7L Water Bottles



14.0416.00 0.7L Water Bottle Service Kit

Includes (1) bottle, tamper-proof bottle cap, washer, supply tube, and barb fitting with washer

14.0332.01 0.7L Water Bottle Supply Tub (6)

004.005.02 Flat Nylon Washer 0.187 ID (10)

023.001.03 1/4 x 10-32 Barb (10)

004.137.00 Flat Urethane Washer 0.750 ID

052.015.00 Tamper-proof Bottle Cap

Replacement Receptacle Assembly

for Self-Contained Water Systems

14.0408.00 Cap assembly


Cascade/Decade/Performer

14.0308.02 Cap assembly


Pre-Cascade Systems

14.0492.00 Cap assembly

Radius w/short stud




Cascade/Decade/
Radius/Performer



Pre-Cascade Systems

Regulator, Relieving, 40 psi

24.0388.02



- For self-contained water system
- Non-adjustable
- 10 - 32 inlet and outlet ports
- Includes 023.004.00 barb
- Includes 004.005.00 washer

A-dec ICX® Waterline Treatment Tablet


90.1064.00 Case, 0.7-liter tablet




- For 0.7-liter bottles
- Case contains (36) boxes with (50) tablets each

Mounting Bracket For Self-Contained Water


14.0327.01




Diaphragmes et joints




Black
24.0137.01 (10)




Black
38.0711.01 (5)




Black
22.0440.02 (10)
22.0440.01 (50)




Black
24.0136.01 (10)




Black
38.0179.02 (10)




Clear
38.0507.01 (10)




Clear
35.0574.01 (5)




Black
38.1787.01 (5)



Black
38.0519.01 (10)




Clear
29.0050.01 (5)




Clear
38.0186.02 (10)
38.0186.01 (25)

NOTE: Number in parentheses
indicates quantity.

Diaphragm and gasket
illustrations are not to scale.



Black
38.0054.02 (10)
38.0054.01 (50)



Black
75.0911.01 (5)

Joint torques

Use A-dec O-Rings


A-dec O-rings are made from specific ethylene propylene (EP), fluorocarbon (Viton) and nitrile (Buna N) compounds to provide the best performance. We recommend the use of A-dec O-rings in your A-dec equipment.

A-dec's silicone lubricant is a high quality silicone-based grease. It is ideal for lubricating A-dec O-rings. Before installing O-rings, always apply a light coating of silicone lubricant. This makes the installation of the O-ring easier and prevents the O-ring from being damaged. A-dec silicone lubricant is available in a package of 6 as Part Number 98.0090.01.

The use of petroleum-based lubricants (jellies, creams, or lotions) cause the O-rings in the syringe to swell and/or deteriorate, resulting in syringe malfunction or failure.

O-ring names

O-ring names, which appear on acknowledgements and invoices, provide the following informative descriptions:



Autoclavable O-Ring Specification Chart Package of 10

PART NUMBER	MATERIAL	I.D.	O.D.	WIDTH	PRICE
034.001.01	EP	.032 in	.108 in	.038 in	
034.002.01	EP	.042 in	.142 in	.050 in	
033.003.01	Viton	.056 in	.176 in	.060 in	
034.003.01	EP	.056 in	.176 in	.060 in	
035.050.01	EP	.070 in	.150 in	.040 in	
035.033.01	EP	.101 in	.241 in	.070 in	
035.030.01	Viton	.176 in	.256 in	.040 in	
035.049.01	EP	.176 in	.256 in	.040 in	
034.107.01	EP	.206 in	.412 in	.103 in	
034.012.01	EP	.364 in	.504 in	.070 in	
034.013.01	EP	.426 in	.566 in	.070 in	
034.014.01	EP	.489 in	.629 in	.070 in	
034.018.02	EP	.739 in	.879 in	.070 in	
034.019.01	EP	.801 in	.941 in	.070 in	

Jointes toriques

Non-Autoclavable O-Ring Specification Chart

PART NUMBER	MATERIAL	I.D.	O.D.	WIDTH	PKG QTY	PRICE
030.001.02	Buna-N	.029 in	.109 in	.040 in	10	
030.002.02	Buna-N	.042 in	.142 in	.050 in	10	
030.003.02	Buna-N	.056 in	.176 in	.060 in	10	
030.004.02	Buna-N	.070 in	.210 in	.070 in	10	
030.005.02	Buna-N	.101 in	.241 in	.070 in	10	
030.008.02	Buna-N	.176 in	.316 in	.070 in	10	
030.009.01	Buna-N	.208 in	.348 in	.070 in	25	
030.011.02	Buna-N	.301 in	.441 in	.070 in	10	
030.012.02	Buna-N	.364 in	.504 in	.070 in	10	
030.013.02	Buna-N	.426 in	.566 in	.070 in	10	
030.014.01	Buna-N	.489 in	.629 in	.070 in	25	
030.015.02	Buna-N	.551 in	.691 in	.070 in	10	
030.016.02	Buna-N	.614 in	.754 in	.070 in	10	
030.019.03	Buna-N	.801 in	.941 in	.070 in	10	
030.021.01	Buna-N	.926 in	1.066 in	.070 in	2	
030.023.02	Buna-N	1.051 in	1.191 in	.070 in	10	
030.026.01	Buna-N	1.239 in	1.379 in	.070 in	10	
030.027.01	Buna-N	1.301 in	1.441 in	.070 in	10	
030.034.02	Buna-N	2.114 in	2.254 in	.070 in	10	
030.035.02	Buna-N	2.239 in	2.379 in	.070 in	10	

Non-Autoclavable O-Ring Specification Chart Con't

PART NUMBER	MATERIAL	I.D.	O.D.	WIDTH	PKG QTY	PRICE
030.224.01	Buna-N	1.734 in	-	.139 in	10	
030.228.02	Buna-N	2.250 in	2.500 in	.125 in	10	
031.130.01	Buna-N	1.612 in	1.818 in	.103 in	10	
034.004.01	EP	.070 in	.210 in	.070 in	10	
035.003.01	Buna-N	.029 in	.109 in	.040 in	10	
035.004.01	Buna-N	.240 in	.306 in	.033 in	10	
035.024.01	Buna-N	.176 in	.256 in	.040 in	10	
035.025.01	Buna-N	.070 in	.150 in	.040 in	10	
035.032.01	Buna-N	.237 in	.317 in	.040 in	10	
035.034.01	Buna-N	.100 in	.176 in	.038 in	10	
035.041.02	Buna-N	.614 in	.768 in	.077 in	10	
035.053.01	Silicone	20 mm	-	2 mm	10	

Caches et couvercles A-dec 500

Accessory Plate

- 77.0147.00 Plate, accessory Surf 4
- 77.0147.01 Plate, accessory White 2



- For A-dec 532 and A-dec 542 Traditional and A-dec 533 Continental delivery systems

Monitor Post Cap

- 28.1501.00 Surf 4
- 28.1501.01 White 2



- Fits the top of the post on an A-dec 561 support-side mounted monitor


Support Link Bearing Cover Cap and Locking Knob

- 77.0374.00 Bearing cover cap, Surf 4
- 77.0374.01 Bearing cover cap, White 2
- 77.0373.00 Knob, link lock, Surf 4
- 77.0373.01 Knob, link lock, White 2



Blank Whip Cover

- 77.0069.00 Surf 4
- 77.0069.01 White 2



- For A-dec 533 Continental delivery systems

Cover, Brake Handle Hole

- 77.0424.00 Surf 4
- 77.0424.01 White 2



- Covers the brake handle hole on A-dec 532 and A-dec 542 Traditional and A-dec 533 Continental delivery systems
- For delivery systems after June 2005

Assistant's Arm Covers

- 77.0021.01 Covers (4), Surf 4
- 77.0021.03 Covers (4), White 2



- Covers the pivot points on an A-dec 551 chair-mounted assistant's arm

Caches et couvercles A-dec 300

Plug, Control head touchpad hole

77.0845.00



- Fits either side of the control head

Plug, Monitor cable hole in light post

75.0150.00



- Plugs light post when no monitor is present

Chair upper structure Hole plug

018.002.02 Gray
018.002.04 Surf 4
018.002.05 Gray 3
018.002.06 White 2
018.002.08 Slate



Light Flexarm end cap

75.2018.00



- Holds additional light bulb
- Replacement bulb 041.709.00
- For A-dec 371 lights prior to serial number F052936

Trim Ring, telescoping assistant arm

77.0877.00



- 2" post trim cover for assistant arm

Indicator Light Cover

018.119.00



- For A-dec Motor/pump cover
- 1" diameter

Vacuum canister blank cover

77.0724.00



- Use on A-dec 361 support center when no vacuum is needed

Umbilical hole cover

62.0348.00



- Use on A-dec 311 Utility cover

End Cap for A-dec 300 Flexarm

77.0784.00



- For all A-dec 300 delivery system flexarms

Caches et couvercles A-dec 300

Cover for Telescoping Assistant's Arm

77.0878.00



- For A-dec 352 telescoping assistant's arm with no light

Cap for 2" post on A-dec 300 Monitor

77.0842.00



- For A-dec 300 monitor mount with no light
- Fits on top of intermediate post

Plug kit for a-dec 300 Support Center

90.1215.00



- For A-dec 361 support center with no light or delivery

End Cap For a-dec 361 water bottle arm

90.1212.00



- For A-dec 361 water bottle arm

Embouts et manchons

1/8" OD In-Line Barb


023.019.03 (10) Barbs

023.019.02 (25) Barbs




5/16" OD In-Line Barb

023.814.01 (5) Barbs



Connector Kit


023.025.00



- (3) Washers, .187 ID
- (1) 10 - 32 plug
- (1) 1/4" barb
- (1) 1/8" barb
- (1) 1/4" tee body
- (1) 1/8" uni-clamp
- (3) 1/4" ID sleeve clamps

1/4" OD In-Line Barb

023.045.04 (10) Barbs



Connector Kit


023.079.00



- (3) Washers, .187 ID
- (3) 1/4" ID sleeve clamps
- (1) 5/16" ID sleeve clamp
- (1) 10-32 Female coupling
- (1) 1/8" uni-clamp
- (1) 1/4" Tee body
- (1) 1/4" x 10-32 Barb
- (1) 1/8" x 10-32 Barb
- (1) 5/16" x 10-32 Barb

1/8" to 1/4" Barb Adapter Kit


023.029.00 Kit



- Includes:
- (5) Washers, .187 ID
- (1) Hex plug, 10 - 32
- (1) Female coupling, 10 - 32
- (2) Barbs, 1/4" x 10 - 32
- (2) Barbs, 1/8" x 10 - 32
- (2) Uni-clamps, 1/8"
- (2) Sleeve clamps, 1/4" ID

Barb/Fittings Service Parts Kit

90.1021.00 Service parts kit



Includes (quantities shown in parentheses):

- (40) 1/8" tubing plug
- (20) 1/4" tubing plug
- (50) 1/8" barb cap
- (40) 1/4" x 10 - 32 barb
- (60) 1/8" uni-clamp
- (20) 1/4" barb cap
- (30) 1/8" x 10 - 32 barb
- (20) Sleeve clamp
- (10) 1/4" tee body
- (30) 1/4" in-line barb
- (30) 1/8" barbed tee
- (10) 1/8" MPT to 10 - 32 adapter
- (60) 1/8" in-line barb
- (50) 10 - 32 hex plug
- (5) 10 - 32 cross body
- (30) 10 - 32 female coupling
- (120) Flat nylon washer
- (5) 10 - 32 to 3/8" compression adapter

300/500 Barb, Fitting Service Kit

90.1268.00




- Includes:
- (4) Barb, 5/16" x 10 - 32
- (2) Barb, 1/8" x 10 - 32
- (2) Barb, 1/4" x 10 - 32
- (4) Flat nylon washers, .187
- (4) Sleeve Clamp, 5/16"
- (2) Uni-clamp, 1/8"
- (1) Female cross, 10 - 32
- (2) Hex plug, 10 - 32

Embouts et manchons

10 - 32 Female Cross Body

021.027.01 (5) Cross bodies



5/16" Water Tubing Barb


023.805.01 Barbs, 10 - 32 (10)

023.804.00 Barb, 1/8" NPT



5/16" x 10 - 32


- Requires washer 004.005.02 (10)



5/16" x 1/8" NPT

1/8" x 1/4-32 Barb


023.028.00 Barb



- Requires washer 004.008.01 (10)

1/8" OD 4-Way Barbed Plastic Tee

023.065.01 (5) Tees



1/8" x 10 - 32 Barb

023.004.03 (10) Barbs

023.004.04 (100) Barbs



- Requires washer 004.005.02 (10)

1/8" x 10 - 32 Restrictor Barb

023.070.00 (1) Barb



.013

- Use in self-contained water system
- Requires washer 004.005.02 (10)


1/8" OD 3-Way Barbed Plastic Tee

023.014.02 3-way (10) Tees



1/8" x 10 - 32 Restrictor Barb

023.036.00 (1) Barb




.005

- Requires washer 004.005.02 (10)

1/4" x 10 - 32 Barb

023.001.03 (10) Barbs

023.001.02 (25) Barbs



- Requires washer 004.005.02 (10)

Embouts et manchons


1/8" Tubing Plug

023.030.02 (10) Plugs



1/4" Tubing Plug

023.057.01 (10) Plugs



1/8" ID Uni-clamp

025.007.01 (10) Uni-clamps

025.007.02 (100) Uni-clamps




- Fits all A-dec 1/8" polyurethane tubing
- Plastic, white

Quick Connector

023.094.01 1/8" (5)

023.095.01 1/4" (5)



- Nickel-plated brass

Sleeve Clamps


025.015.00 1/4" Tubing clamp (1)

025.015.01 1/4" Tubing clamps (10)

025.015.02 1/4" Tubing clamps (100)

025.109.00 5/16" Tubing clamp (1)

025.109.01 5/16" Tubing clamps (10)



025.015.XX

- Plastic

025.109.00
025.109.01

10 - 32 Hex Plug

021.016.04 (10) Plugs




- Nickel-plated brass
- Requires washer 004.005.01

Barb Cap

023.031.01 1/8" OD (10)

022.009.01 1/4" OD (10)



- Requires 1/8" barb cap uni-clamp 025.007.01

10 - 32 Female Coupling

021.034.01 (10) Couplings



- Nickel-plated brass

.187 ID Washer


004.005.01 (100) Washers

004.005.02 (10) Washers



- Nylon
- Use with 10 - 32 threaded fittings

Attaches rapides




Quick Disconnect Fittings for Self-contained Water

13.0375.00	Male	1/4" OD with 1/4" barb	
026.125.00	Female	1/4" ID with 1/4" & 1/8" barb	
90.0370.00	Kit*		

A-dec 500 Water QD Kit

90.1097.00 Water QD kit




Quick Disconnect Fittings for City Water

13.0381.00	Male	1/4" OD with 1/4" barb	
026.131.00	Female - Cascade	1/4" ID with 1/4" barb	
026.157.00	Female - Performer	1/4" ID with 1/8" barb	
026.134.00	Kit* - Cascade		
76.2025.00	Kit* - Performer		

A-dec 500 Air QD Kit

90.1098.00 Air QD kit




Quick Disconnect Fittings for Air

13.0384.00	Male	1/4" OD with 1/4" poly	
026.132.00	Female	1/4" ID with 1/4" barb	
026.130.00	Kit* - Cascade		
76.2030.00	Kit* - Performer		

A-dec 300 Air and Water QD Kit

90.1166.00 Water QD kit



*Kits include male and female connectors and required hardware

NOTE: Male and female connectors only connect with items in the same table, they cannot be mixed and matched

Profilées

Boîtes au sol 311, 411, 511

		• Cast frame • Does not include umbilical or utilities • 17" x 16.5" (432 mm x 419 mm) • With outer convolute umbilical cover	
PRODUCT DESCRIPTION	PART NUMBER SURF	PART NUMBER WHITE	PRICE
Floor box and cover	41.1506.00	41.1506.01	
Floor box frame only	41.1467.00	41.1467.01	
Cover only	41.1464.01	41.1464.02	
Master toggle kit	90.1088.00	90.1088.01	
Air filter/regulator	24.0475.00		
Water filter/regulator	24.0474.00		
18" Outer Convolute umbilical cover	024.203.01	024.228.01	

Cascade Contoured Floor Box

	• Cast frame • Surf 4 cover • Does not include umbilical or utilities • 14.5" x 17" (368 mm x 432 mm)	
PRODUCT DESCRIPTION	PART NUMBER	PRICE
Cascade contoured floor box frame and cover	41.1457.00	
Floor box frame only	41.0417.00	
Cover only, Surf 4	41.0416.00	

A-dec 300 avant 2015

	• White cover • Does not include umbilical or utilities • Frame and cover kit also contains all low chair covers	
PRODUCT DESCRIPTION	PART NUMBER	PRICE
Floor box frame only	62.0658.00	
Cover only	62.0655.00	

Mounting Kit


For use with

- Stainless steel floor boxes
- Cascade contoured floor box

41.1454.01 White

Arm Bearing Kit

90.1123.00 Bearing Stack
 016.108.00 Race Thrust Bearing (1)
 016.044.00 Needle Thrust Bearing (1)
 048.107.00 Lubriplate



- Kit includes (2) Race Bearings

Inoxidables

A-dec 500 Stainless Floor Boxes



- Includes 7' (2,134 mm) umbilical
- Use with non-chair mounted delivery system
- 15" x 7" x 17.5" (381 mm x 178 mm x 445 mm)
- Includes mounting hardware
- Includes stainless steel frame and cover

PRODUCT DESCRIPTION	PART NUMBER SURF	PART NUMBER WHITE	PRICE
Remote stainless steel floor box with cover for 541/545	-	41.1490.01	
Remote stainless steel floor box with cover for 542	-	41.1491.01	
Air filter/regulator	24.0475.00		
Water filter/regulator	24.0474.00		

A-dec 300 Remote Floor Box



- Includes 3' (914 mm) umbilical
- Hole plugs
- Mounting hardware
- 5.5" x 14.5" x 17" (140 mm x 368 mm x 432 mm)

PRODUCT DESCRIPTION	PART NUMBER	PRICE
Floor box frame only	41.1507.00	
Cover only, Gray 3	41.0413.01	
Remote stainless steel floor box with cover for 300	77.1019.00	

Large Stainless Steel Floor Box



- Does not include umbilical or utilities
- 5.5" x 14.5" x 17" (140 mm x 368 mm x 432 mm)

PRODUCT DESCRIPTION	PART NUMBER	PRICE
Large floor box frame and cover, Surf 4	41.1456.00	
Large floor box frame and cover, Gray 3	41.1456.01	
Floor box frame only	41.0414.00	
Cover only, Surf 4	41.0413.00	
Cover only, Gray 3	41.0413.01	


Medium Stainless Steel Floor Box



- Surf 4 cover
- Does not include umbilical or utilities
- 5.5" x 9" x 14" (140 mm x 229 mm x 356 mm)

PRODUCT DESCRIPTION	PART NUMBER	PRICE
Medium stainless steel floor box frame, cover and mounting kit	41.1455.00	
Floor box frame with plugs	41.0408.00	
Cover only, Surf 4	41.0407.00	

Performer Stainless Steel Floor Box



- Surf 4 cover
- Does not include umbilical or utilities
- 5.5" x 13" x 13" (140 mm x 330 mm x 330 mm)

PRODUCT DESCRIPTION	PART NUMBER	PRICE
Floor box frame only	41.1317.00	
Cover only, Surf 4	41.0016.02	
Indicator light assembly	041.582.00	

Bouchons et tuyau ombilical


1/2" Hole Plug

018.002.02 Gray
018.002.04 Surf 4
018.002.05 Gray 3
018.002.06 White 2
018.002.08 Slate




1-1/4" Hole Plug

018.010.02 Surf 4
018.010.03 Gray 3




Retaining Clips

41.0074.00 1 3/4"
41.0411.00 2 3/8"
47.1259.00 2 3/8" for contour floor box




2" Hole Plugs

018.026.02 Surf 4, round
018.026.03 Gray 3, round
018.026.04 White 2, round
018.026.05 Slate, round
41.1299.00 Surf 4




Umbilical Restraints

41.1143.00 Straight restraint
41.1144.00 U-shaped restraint



Umbilical Connectors


41.1199.00 1" Assembly
75.0948.00 1 1/2" Assembly
41.0062.01 1 1/4" Connector
41.1164.00 1 3/4" Connector
41.1385.00 2" Connector



• Surf 4 color

Chair Cord Bushing

41.0192.02 Surf 4
41.0192.03 Gray 3



Screw


003.011.01 (5)



- #10 sheet metal pan head, 11/4", Zinc
- Mounting hardware for floor boxes

Mounting Hardware

014.007.00 Toggle wall anchor
014.006.00 Threaded insert



1/4-20 Wall anchor toggle bolt



3/16 x 2" Threaded Wall Insert

Séparateurs d'eau et filtres

Moisture Separator, Manual

026.186.00 Moisture separator


97.0405.00 Bowl assembly



- 1/8" MPT ports
- Bowl assy 97.0405.00
- After December 2005

Moisture Separator Bowl with Seal

97.0290.00 Bowl with seal only




- For moisture separator 026.066.00 and automatic moisture separator 90.1027.03 before December 2005
- Uses filter 97.0280.02 (6)

Filter, Moisture Separator, Automatic

90.1027.03 Universal kit

97.0406.00 Bowl assembly


97.0408.01 Filter (3)



- 1/8" MPT ports
- After December 2005

Air Filter/Regulator/Moisture Separator

026.021.01



- Air pressure regulator and filter combination
- Manually drained
- Includes gauge, fitting, mounting bracket, and hardware

Filter For Moisture Separators

97.0280.02 (6)

97.0408.01 Filter (3)

97.0401.01



- 5-micron filtration (not a bacterial filter)
- For moisture separators 026.066.00, 41.0217.00, 47.1147.03 use 97.0280.02
- For air filter/regulator/moisture separator 026.021.01 use 97.0280.02 (before 2003) and 97.0401.01 (beginning 2003)

A-dec 500 Moisture Separator


41.1477.00 Moisture sep. assy



Air Filter/Regulator Service Kit

90.0030.00 Before 2003

90.1066.00 Beginning 2003




- Kit to repair the air filter regulator
- Includes diaphragm, poppet, spring, filter element and short valve core

InLine Filter

19.0100.00 Assembly

19.0220.01 Filters (10)




- For air or water; 1/4" barb fittings
- Includes 5/8" x 1/16," 25 - 30 micron filter (non-bacterial)

Filtres/régulateurs

Universal Air and Water Filter/Regulator Kit


24.0471.00 Kit



- Includes gauges, adjustable regulators, non-bacterial filters, and pre-regulators
- Works for both Cascade and Performer
- Replaces 24.0363.00 and 24.0363.04

Air Regulator Assembly


24.0128.01



- Works on any A-dec air regulator manifold
- Includes 9-hole gasket and mounting hardware

Universal Air Filter/Regulator


24.0469.00 Filter/regulator



- Includes gauge, adjustable regulator, non-bacterial filter, and pre-regulator
- Works for both Cascade and Performer
- Replaces 24.0363.00 and 24.0363.04

Water Regulator/Shutoff with Spacer


24.0138.01



- 2:1 ratio of pilot air to water outlet pressure
- Includes water regulator with spacer, 9-hole gasket, and mounting hardware

Universal Water Filter/Regulator


24.0470.00 Filter/regulator



- Includes gauge, water regulator, pre-regulator, non-bacterial filter, and 3/8" and 1/4" adapter, 1/4" QD

Water Regulator/Shutoff

24.0129.01



- 1:1 ratio of pilot air to water outlet pressure
- Includes regulator and gasket

Filtres, pré-régulateurs, manomètres, valve d'arrêt

Filter

24.0234.01 (6)



- Filters air or water
- 70 – 90 micron filtration (not a bacterial filter)

Gauge, Panel Mount Kit

026.118.00 Kit




- Panel mount
 - 0 – 100 psi
 - 1 1/2" diameter
 - 1/8" MPT x 10 – 32 threads
- Includes bracket and hardware

Manual Shutoff Valve

026.062.00 1/2" compression

026.158.00 1/2" female pipe thread




026.158.00

- Inlet 1/2" Compression or 1/2" Female pipe thread
- Outlet 3/8" compression
- Nickel-plated brass


Air or Water Regulator, Non-Relieving

24.0251.01



Pre-Regulator, Relieving


24.0182.02



- Controls air and water pressure in most A-dec systems built since 1978
- Female 10 – 32 threads
- Includes 006.009.00 nut

Regulator Assembly Pre-Set 40 psi

24.0382.00



- Preset to 40 psi in the factory, not adjustable
- 1/8" MPT or 10 – 32 inlet port, 10 – 32 outlet
- Replaces 24.0118.00

Clés hexagonales, test anti-rétraction, manomètres

Panel Mount Gauge


026.118.00



- 0 - 100 psi
- 1 1/2" diameter
- Used to check air/water pressure valves at various system locations
- 1/8" MPT x 10 - 32 threads

Water Retraction Tester


98.0051.00



- Use to determine whether or not a handpiece control system retracts water
- Includes instructions

Handpiece Pressure Gauge


50.0271.00



- 0 - 60 psi
- Fits all handpiece couplers except Borden 3-hole

Hex Wrench Set

009.018.00




- Includes hex wrenches required in servicing A-dec equipment
- Plastic holder

Cable Ties

025.036.01 (10) Small ties

025.002.01 (10) Medium ties

025.003.01 (10) Large ties




- Three sizes:
Small: 3 3/4" (9.5 cm)
Medium: 6 3/4" (17 cm)
Large: 14" (35 cm)

Ball End, Hex Driver

009.037.00 1/8," Small handle

009.005.00 3/32," Std handle


008.023.00 5/64," Std handle




Clips and Cable Tie Mount

025.017.01 Cord clip (10)

025.019.01 Mount (10)



Adhesive-backed cord clip



Adhesive-backed cable tie mount

Pac 1 Toolkit

45.0438.00


Track Light Assembly Tool

28.1832.00

Peinture de retouche et silicone

Syringe Tool

98.0010.00



- Use to remove/replace valve cores on Uni-Flo and Tri-Flo Syringes
- Use 3/32" to adjust syringe flow on Century II
- Use 5/32" to remove/replace the syringe tip retaining nut

Valve Test Syringe

98.0050.01 (4)



- For making quick tests of pilot-operated valves
- Applies a static pressure of 5 to 75 psi

A-dec Silicone Lubricant

98.0090.01 0.05oz blister packs (6)

98.1120.00 Re-sealable tube 5.3oz



Touch-Up Paint

90.0351.00 Surf 4

90.0351.01 Surf 6

90.0351.02 White 2

90.0351.03 Gray 3



- 1/2 fluid ounce
- Includes brush



Loctite®

060.001.00 Red #271

060.002.00 Blue #242



Service Tip


A-dec's silicone lubricant is a high quality silicone-based grease. It is ideal for lubricating A-dec O-rings. Before installing O-rings always apply a light coating of silicone lubricant. This makes the installation of the O-ring easier and prevents damage to the O-ring. A-dec silicone lubricant is available in a package of 6 as Part Number 98.0090.01.

Caution

L'utilisation de lubrifiants à base de pétrole provoque une détérioration précoce des joints toriques.

Codification des tubulures 1/8^e

Tubing — Color-Coded Polyurethane, 1/8"



1/8" OD Polyurethane Tubing

Durometer: 90A; working pressure: 150 psi maximum; suggested fittings: 1/8" barbs with sleeves

TUBING COLOR	TRACER COLOR	TRACER MARKING	10' (3 M)	50' (15 M)
Yellow	None	None	036.003.03	036.003.04
Yellow	Red	Stripe	036.013.03	036.013.04
Yellow	Red	Long Dash	036.009.04	036.009.05
Yellow	Green	Long Dash	036.044.03	036.044.04
Orange	Black	Long Dash	036.010.03	036.010.04
Brown	White	Long Dash	036.014.02	-
Green	White	Long Dash	036.006.03	036.006.04
Green	White	Short Dash	036.018.03	036.018.04
Blue	None	None	036.004.03	036.004.04
Red	None	None	036.005.03	036.005.04
Purple	None	None	036.011.03	-
Black	None	None	036.020.03	-
Clear	None	None	024.015.04	024.015.05
White	None	None	036.019.03	036.019.04


1/8" OD Vinyl Tubing

TUBING COLOR	TRACER COLOR	TRACER MARKING	10' (3 M)	50' (15 M)
Dark Surf	None	None	036.043.01	036.043.02

NOTE: Automatic shutoff must be installed at manual shutoff for water applications.

Codification des tubulures suite

Tubing — Color-Coded Polyurethane 1/4," 5/16," 3/8"



1/4" OD Polyurethane Tubing

Durometer: 90A; Working pressure: 100 psi maximum; Suggested fittings: 1/4" barbs with sleeves and 1/4" poly-flo

TUBING COLOR	TRACER COLOR	TRACER MARKING	10' (3 M)	50' (15 M)
Yellow	Black	Long dash	036.051.03	036.051.04
Yellow	None	None	036.103.03	036.103.04
Orange	Black	Long dash	036.052.03	036.052.04
Blue	None	None	036.053.03	036.053.04
Blue	White	Stripe	036.054.03	036.054.04
Clear	None	None	036.066.03	036.066.04
Grey	None	None	024.090.04	024.090.05

5/16" OD Polyurethane Tubing

Durometer: 90A; Working pressure: 100 psi maximum; Suggested fittings: 5/16" barbs with sleeves

TUBING COLOR	TRACER COLOR	TRACER MARKING	10' (3 M)	50' (15 M)
Yellow	None	None	036.114.01	036.114.03
Blue	None	None	036.116.01	036.116.03

3/8" OD Polyurethane Tubing




Durometer: 90A; Working pressure: 100 psi maximum; Suggested fittings: 3/8" poly-flo

TUBING COLOR	TRACER COLOR	TRACER MARKING	10' (3 M)	50' (15 M)
Yellow	Black	Stripe	036.032.02	036.032.04
Yellow	None	None	036.031.02	036.031.04
Blue	White	Long dash	036.033.02	036.033.04
Grey	None	None	036.034.02	036.034.04

Cordons en silicone

Silicone Handpiece Tubing

- Quick-connect style tubing
- 61" (1,549 mm) for use with A-dec 532/533/542 and A-dec 332/333/334/335
- 85" (2,159 mm) for use with A-dec 541

	PRODUCT DESCRIPTION	TUBING LENGTH	PART NUMBER DARK SURF	PART NUMBER GRAY 3	PRICE
	3-hole tubing For use with side exhausting Borden-style tubing	61" (1,549 mm) 85" (2,159 mm)	98.1063.00 98.1068.00	98.1086.00 98.1109.00	
	4-hole tubing	61" (1,549 mm) 85" (2,159 mm)	98.1064.00 98.1068.00	98.1087.00 98.1106.00	
	6-pin fiber-optic	61" (1,549 mm) 85" (2,159 mm)	98.1066.00 98.1070.00	98.1089.00 98.1108.00	

Silicone Handpiece Tubing QD Barbs

026.174.01 QD, 1/8" barb (5)

026.175.01 QD, 3/16" barb (5)



1/8"

- For A-dec 300/500 silicone handpiece tubing
- Includes O-ring



3/16"

Cordons en vinyl


Vinyl Handpiece Tubing

	PRODUCT DESCRIPTION	TUBING LENGTH	PART NUMBER DARK SURF	PART NUMBER GRAY 3	PRICE
	3-hole tubing For use with side exhausting Borden-style tubing	84" (2,134 mm)	98.0882.00	98.1093.00	
	4-hole tubing	84" (2,134 mm) 144" (3,657 mm)	98.0879.00 98.0879.12	98.1094.00 98.1113.00	
	6-pin fiber-optic	84" (2,134 mm) 144" (3,657 mm)	98.0885.00 98.0885.12	98.1096.00 98.1112.00	

Vinyl Handpiece Tubing QD Barbs

026.161.00 QD, 1/8" barb


023.088.00 QD, 1/4" barb



- For attaching vinyl tubing to the A-dec 300/500 control block
- Includes O-ring

Cordons en vinyl, raccords et ampoules

Midwest Handpiece Tubing (straight dark surf)

	PRE-CUT LENGTH	PER FOOT	10' (3 M)	50' (15 M)	100' (30 M)
	3-Hole	024.147.00	-	024.147.01	-
	4-Hole	-	024.148.02	024.148.01	024.148.04

Handpiece Tubing Nuts

- 98.0180.00 Borden handpiece nut
- 98.0190.00 Midwest handpiece nut
- 98.1119.00 Fiber-optic, 5-hole handpiece nut
- 98.1048.00 Fiber-optic, 6-pin handpiece nut




Nut

Terminal Nut Assembly

- 98.0450.00 Midwest handpiece, Terminal Nut Assembly

A-dec 500 et 300 Aspiration grand volume

52" (130 cm) tuyau utilisé pour les 351 & 551 instrumentation montée sur fauteuil
63" (160 cm) tuyau pour les 352, 353, 541 & 545





REF	DESCRIPTION	PART NUMBER	PRICE
Standard (11mm) HVE Valve (for 1/2 inch or 16mm tubing)			
1	Standard (11mm) HVE Valve Assembly	11.1335.00	
2	Tip o-ring (10)	034.013.01	
3	HVE Rotor Assembly	11.1074.00	
4	Rotor o-ring (10)	034.014.01	
5	Spring Clip Spring (5)	11.0998.01	

1/2 inch tubing

REF	DESCRIPTION	SURF	WHITE	PRICE
6	Standard (11mm) HVE Assembly with 1/2 in. tubing, valve, and (2) tailpieces			
	with 52" (1321mm) tubing	11.1286.00	11.1286.01	
	with 63" (1600mm) tubing	11.1296.00	11.1296.01	
7	1/2 inch Tubing with (2) Standard Tailpieces			
	with 52" (1321mm) tubing	11.1292.00	11.1292.01	
	with 63" (1600mm) tubing	12.1237.00	12.1237.01	
	Tailpieces			
8	Tailpiece for Valve End	11.1027.00	11.1027.01	
9	Tailpiece for Vacuum End	11.1300.01	11.1300.03	


16mm tubing

REF	DESCRIPTION	SURF	WHITE	PRICE
10	Standard (11mm) HVE Assembly with 16mm tubing, and (2) tailpieces			
	with 52" (1321mm) tubing	11.1288.00	11.1288.01	
	with 63" (1600mm) tubing	11.1305.00	11.1305.01	
11	16mm Tubing with (1) 11mm Tailpiece and (1) 16mm Tailpiece			
	with 52" (1321mm) tubing	11.1294.00	11.1294.01	
	with 63" (1600mm) tubing	12.1239.00	12.1239.01	
	Tailpieces			
12	11mm Tailpiece for Valve End	11.1272.00		
13	16mm Tailpiece for Vacuum End	12.1121.00	11.1301.03	




A-dec 500 et 300 Pompe à salive

52" (130 cm) tuyau utilisé pour les 351 & 551 instrumentation montée sur fauteuil
63" (160 cm) tuyau pour les 352, 353, 541 & 545



REF	DESCRIPTION	PART NUMBER	PRICE
Standard SE Valve (for standard or 10mm tubing)			
1	Standard SE Valve Assembly	12.1258.00	
2	Tip o-ring (10)	034.107.01	
3	SE Rotor Assembly	12.1093.00	
4	Rotor o-ring (10)	034.012.01	
5	Screen (10)	11.1235.01	




10mm tubing


REF	DESCRIPTION	SURF	WHITE	PRICE
10	Standard SE Assembly with 10mm tubing, Valve, and (2) tailpieces			
	with 52" (1321mm) tubing	12.1228.00	12.1228.01	
	with 63" (1600mm) tubing	12.1236.00	12.1236.01	
11	10mm Tubing with (2) Tailpieces			
	with 52" (1321mm) tubing	12.1232.00	12.1232.01	
	with 63" (1600mm) tubing	11.1307.00	11.1307.01	
Tailpieces				
12	Tailpiece for Valve End	11.1283.01	11.1283.03	
13	Tailpiece for Vacuum End	11.1302.01	11.1302.03	

Cascade et Performer Aspiration grosse canule

72" (180 cm) pour Performer
84" (210 cm) pour les Cascade
144" (360 cm) pour les carts Cascade




REF	DESCRIPTION	PART NUMBER	PRICE
	Standard (11mm) HVE Valve (for 1/2 inch or 16mm tubing)		
1	Standard (11mm) HVE Valve Assembly	11.1335.00	
2	Tip o-ring (10)	034.013.01	
3	HVE Rotor Assembly	11.1074.00	
4	Rotor o-ring (10)	034.014.01	
5	Spring Clip Spring (5)	11.0998.01	



HVE with 1/2 inch tubing


REF	DESCRIPTION	PART NUMBER	PRICE
6	HVE Assembly with 1/2 inch tubing, valve, and (2) tailpieces		
	Standard (11mm) HVE		
	with 84" (2434mm) tubing	11.1025.02	
	with 144" (3658mm) tubing	11.1025.12	
	Long Tip Holder HVE		
	with 84" (2434mm) tubing	11.1178.00	
	with 144" (3658mm) tubing	11.1014.12	
	Performer HVE for AVS		
	with 84" (2434mm) tubing	11.1127.01	
	Tailpieces		
7	Tailpiece for Valve End		
	for Standard (11mm) HVE or Long Tip Holder HVE	11.1027.00	
	for Performer HVE for AVS	11.1131.00	
8	Tailpiece for Vacuum End		
		11.1027.00	

REF	DESCRIPTION	PART NUMBER	PRICE
	Standard SE Valve (for standard or 10mm tubing)		
9	Standard SE Valve Assembly	12.1258.00	
10	Tip o-ring (10)	034.107.01	
11	SE Rotor Assembly	12.1093.00	
12	Rotor o-ring (10)	034.012.01	
13	Screen (10)	11.1235.01	




SE with Standard Tubing

REF	DESCRIPTION	PART NUMBER	PRICE
14	SE Assembly with Standard tubing, Valve, and (1) tailpiece		
	Standard SE		
	with 84" (2434mm) tubing	12.0910.06	
	with 144" (3658mm) tubing	12.0910.12	
	Performer SE for AVS		
	with 72" (1829mm) tubing	12.1246.00	
	Tailpieces		
15	Tailpiece for Valve End		
	for Standard SE	12.1088.00	
	for Performer SE for AVS	12.1087.00	




A-dec 500, 300 et cascade valves d'aspiration de 15 mm

52" (130 cm) tuyau utilisé pour les 351 & 551 instrumentation montée sur fauteuil
63" (160 cm) tuyau pour les 352, 353, 541 & 545



REF	DESCRIPTION	PART NUMBER	PRICE
Large Bore (15mm) HVE Valve (for 16mm Tubing)			
1	Large Bore (15mm) HVE Valve Assembly	12.1257.00	
2	HVE Rotor Assembly	12.1116.00	
3	Rotor o-ring (10)	034.019.01	
4	Tail end o-ring (10)	034.018.02	
5	Spring Clip Spring (5)	12.1109.01	



16mm tubing


REF	DESCRIPTION	SURF	WHITE	PRICE
6	Large Bore (15mm) HVE Assembly with 16mm tubing, valve, and (2) Tailpieces			
	with 52" (1321mm) tubing	11.1287.00	11.1287.01	
	with 63" (1600mm) tubing	11.1297.00	11.1315.00	
	with 84" (2134mm) tubing	12.1132.00		
7	16mm Tubing with (2) 16mm Tailpieces			
	with 52" (1321mm) tubing	11.1293.00	11.1293.01	
	with 63" (1600mm) tubing	12.1238.00	12.1238.01	
	16mm Tailpieces			
8	Tailpiece for Valve End	12.1121.00	11.1301.03	
9	Tailpiece for Vacuum End	12.1121.00	11.1301.03	


Tuyaux d'aspiration, pédale au sol, boîte au sol et crachoir

Tuyau pompe à salive, Vinyl		LENGTH	DARK SURF	GRAY 3	PRICE
	10mm ID (Int'l elect holders)	10' (3m)	024.197.01	024.219.01	
		25' (7.6m)	024.197.02	024.219.02	

Grosse canule, Vinyl		LENGTH	DARK SURF	GRAY 3	PRICE
	5/8" (16mm) ID (for 11mm & 15mm Assemblies)	10' (3m)	024.190.01	024.222.01	
		50' (15m)	024.190.02	024.222.02	

Vacuum and Umbilical Tubing		LENGTH	DARK SURF	GRAY 3	PRICE
25mm ID Vacuum Tubing (for A-dec 200 and 300)		8' (2.5m)	-	024.209.02	
		12' (4m)	-	024.209.01	
20mm ID Cuspidor Drain Tubing (for A-dec 200 and 300)		7' (2m)	-	024.210.02	
		12' (4m)	-	024.210.01	
25mm ID Vacuum Tubing (for A-dec 500)		9' (2.7m)	-	024.244.02	
		25' (8m)	-	024.244.03	
20mm ID Drain Tubing (for A-dec 500)		10' (3m)	-	024.198.01	
		25' (8m)	-	024.198.02	
2 1/2" (63.5mm) ID Umbilical Tubing		10' (3m)	024.203.02	024.228.02	
1/2" (12.7mm) ID Cuspidor Drain Tubing (for Performer)		10' (3m)	024.162.01	-	
		50' (15m)	024.162.03	-	
5/8" (15.9mm) ID Cuspidor Drain Tubing		10' (3m)	024.177.01	-	
		50' (15m)	024.177.03	-	
3/4" (19mm) ID Cuspidor Drain Tubing		10' (3m)	024.176.06	-	
		50' (15m)	024.176.08	-	


Foot Control Tubing	LENGTH	NO SHEATH	DARK SURF	GRAY 3	PRICE
 Tubing Without Sheath	10' (3m)	024.151.01	-	-	
	50' (15m)	024.151.03	-	-	
 Foot Control III Tubing	10' (3m)	-	024.179.01	-	
	50' (15m)	-	024.179.03	-	
 A-dec 500 Foot Control Tubing	10' (3m)	-	024.202.01	024.229.01	
	50' (15m)	-	024.202.02	024.229.02	
 A-dec 300 Foot Control Tubing	10' (3m)	-	-	024.214.01	
	50' (15m)	-	-	024.214.02	

Convolute Tubing for Vacuum or Umbilical				
 Inside Diameter	LENGTH	DARK SURF	GRAY 3	PRICE
1/2" (12.7mm)	10' (3m)	024.193.01	024.220.01	
	50' (15m)	024.193.03	-	
5/8" (15.9mm)	10' (3m)	024.152.01	-	
	50' (15m)	024.152.03	-	
3/4" (19mm)	10' (3m)	024.192.01	-	
	50' (15m)	024.192.03	-	
1" (25.4mm)	10' (3m)	024.156.09	-	
	50' (15m)	024.156.11	-	
1 1/4" (31mm)	10' (3m)	024.184.01	-	
	50' (15m)	024.184.03	-	
1 1/2" (38mm)	10' (3m)	024.166.01	-	
	50' (15m)	024.166.03	-	
1 3/4" (45mm)	10' (3m)	024.153.18	024.232.01	
	50' (15m)	024.153.20	-	
2" (50mm)	10' (3m)	024.191.05	-	
	50' (15m)	024.191.07	-	

Momentanées

2-Way Push-Button

33.0077.01




- Use with air or water
- Panel mount 1/2" hole (1/4" max. panel thickness)
- 1/8" barbs
- Metal button

006.069.00

006.009.00

2-Way Toggle

33.0008.02




006.009.00

- Use on cuspidor cup fill spout
- Panel mount 1/2" hole (1/4" max. panel thickness)
- 1/8" barb, 10 - 32 female coupling
- Surf 4 plastic toggle

2-Way Push-Button

33.0076.01




006.069.00

006.009.00

- Use for air or water
- Panel mount 1/2" hole (1/4" max. panel thickness)
- 10 - 32 rear ported
- Metal button

2-Way Toggle

33.0009.04




006.069.00

006.009.00

- Use for air or water
- Panel mount 1/2" hole (1/4" max. panel thickness)
- 1/8" barbs
- Surf 4 plastic toggle

2-Way Toggle

33.0034.02



006.069.00

006.009.00

- Use for air or water
- Panel mount 1/2" hole (1/4" max. panel thickness)
- 10 - 32 ported
- Surf 4 plastic toggle

Flush Toggle Valve

33.0168.00 Flush toggle Surf 4


33.0168.01 Flush toggle Gray 3



- For A-dec 500 and 300
- 4-way
- 1/8" barbs

2-Way Push-Button


33.0108.02



- Use on Decade cuspidor bowl rinse and cup fill
- Panel mount
- 10 - 32 ports
- Metal button

2-Way Toggle

42.0683.04



006.069.00

006.009.00

- Use for air or water
- Panel mount 1/2" hole (1/4" max. panel thickness)
- 10 - 32 rear ported
- Surf 4 plastic toggle


Manuel ON/OFF, réglage spray

NOTE

Use two-way (non-exhausting) control valves for air, water and cuspidor bowl rinse and cup fill applications. Use three-way (exhausting) control valves for the master toggle, air coolant and water coolant signal applications.

3-Way Toggle


33.0048.05 Surf 4 toggle



- Use with air signal valve
- Panel mount 1/2" hole (1/4" max. panel thickness)
- 1/8" barbs
- Surf 4 plastic toggle

3-Way Sub-Miniature Toggle

026.123.00




- Use with air signal valve
- Panel mount 1/4" hole (1/8" max. panel thickness)
- 1/8" barbs
- Non-adjusting metal toggle

2-Way Toggle Valve Kit

3-Way Toggle Valve Kit

33.0096.02 2-way


33.0080.03 3-way



- Panel mount 1/2" hole (1/4" max. panel thickness)
- 10 - 32 ports
- Surf 4 plastic toggle
- 2-way, air or water
- 3-way, air

2-Way Toggle


12.1070.01



- Use for Performer
- 1/8" MPT ports
- Surf 4 toggle

3-Way Exhausting Toggle


12.1071.02



- Use for Performer
- Surf 4 toggle

Water Coolant

38.0098.00



- For Century I control block
- Manifold mounting

Master Toggle Valves

33.0131.02 For Cascade/Radius Surf 4

33.0173.00 For A-dec 500 and 300 Gray 3

90.1189.00 For A-dec 300 with no delivery

33.0048.05 For A-dec 500 Surf 4




- On/Off toggles

Robinet de réglage

Needle Valve, Key


12.0920.00



- Use for air applications
- For Cascade cuspidor-timed circuits
- Manifold mounting
- 10 – 32 male inlet
- 10 – 32 female outlet

Needle Valve, Lo-Flo


12.0954.01



- Use for air applications
- Panel mount 21/64" hole
(1/8" max. panel thickness)

Flow Control Needle Valve Fine Adjustment

33.0084.01




- Use for coolant water
- Panel mount 1/2" hole
(1/4" max panel thickness)
- 10 – 32 rear ports
- Plastic black stem

Control Valve Assemblies

13.0361.00 Plastic stem


13.0350.00 Brass stem



- Use for air or water applications
- Panel mount 1/2" hole
(1/4" max. panel thickness)
- 10 – 32 ports

Flow Control Needle Valve Coarse Adjustment

33.0085.02




- Use for air or water
- Panel mount 1/2" hole
(1/4" max panel thickness)
- 10 – 32 rear ports
- White plastic stem

Relais

Time Bowl Rinse Relay

12.0913.00




- For Cascade cuspidor
- 3/8" female Poly-Flo
- 1/8" barb
- 3/8" x 1/4" barb

Cup Fill Relay and Cup Fill Timer Stabilizer Kit

12.0934.00 Cup fill relay

90.0456.00 Stabilizer kit



12.0934.00




90.0456.00

- For Cascade cup fill
- 1/8" and 1/4" barb

NOTE: Kit replaces the cup fill timer and improves flow adjustment.

Time Delay Relay


42.0534.00



- For filter regulator manifolds and Excellence cuspidor
- 1/8" and 1/4" barb
- 9-hole gasket, 24.0137.01 (10)

Relay Valve with Holdback

34.0043.00



- For tooth dryer
- Non-Cascade
- 1/8" barbs (1)
- 1/4" barbs (2)

Service Tip


Air bubbles in the water or sputtering may be caused by:

- Loose syringe handle
- Defective O-ring on the connector tube at the syringe terminal
- Leaking diaphragm in the water regulator
- Leaking diaphragm in any water relay valve

Coupure d'eau

Air-Actuated Water Shutoff

34.0031.00




- Inlet 3/8" male compression fitting
- 1/8" barb

Air-Actuated Water Shutoff

34.0032.00 1/4" Poly-flo


34.0033.00 3/8" Poly-flo



- Inlet 3/8" male compression fitting
- 1/4" poly sleeve with 1/4" poly adapter

Air-Actuated Water Shutoff

34.0036.00




- Inlet and outlet 3/8" female compression fittings
- For Preference dental furniture sink

Air-Actuated Air Shutoff

34.0037.01 1/4" Poly-flo

34.0037.00 3/8" Poly-flo




- Inlet 3/8" male compression fitting

Manual Shutoff

026.062.00 1/2" Compression

026.158.00 1/2" Female pipe thread



026.158.00 shown

- Inlet 1/2" compression or 1/2" female pipe thread
- Outlet 3/8" compression
- Nickel plated brass


Service Tip

Restricted water flow could be the result of a clogged screen in the manual shutoff valve. To check:

1. Turn off the master On/Off valve toggle.
2. Bleed off system pressure.
3. Turn OFF manual shutoff valve.
4. Remove water filter/regulator.
5. Check screen.

3-Way Selector


33.0104.05



- Air applications only
- 2 position, 3-way
- 1/8" barbs
- Surf 4 plastic toggle

Ported 3-Way

82.0305.04




- Air or water applications
- 10 - 32 ports (3)
- Panel mount
- Surf 4 plastic toggle

Valves de sélection

2- or 3-Way Micro

33.0134.00 2-way (brass)

33.0138.00 3-way (SST)




- 2-way, (brass ball) for Cascade bowl rinse and Excellence
- 3-way, stainless steel ball (SST), for Cascade cup fill and Excellence (replaces 29.0185.00)

2- or 3-Way Micro

38.0086.00 2-way

29.0146.00 3-way



- 2-way, balanced for vacuum cuspidor
- 3-way, reverse

In-Line Check


026.073.02



- 1/8" tubing barbs
- 15 psi minimum
- Arrow indicates direction of flow

Shuttle

026.074.01



- 1/8" tubing barbs

Tee Check

27.0005.03



- For Century I and II handpiece flush
- 1/8" barbs
- Arrow indicates direction of flow

2-Way Air Bleed

33.0132.00



- For Cascade unitized holder
- Use for air bleed
- Snap fit mounting
- 1/8" barb

2-Way Air Bleed

33.0025.00



- For Excellence holder and Cascade individual holder
- Use for air bleed
- Snap fit mounting
- 1/8" barb

2-Way Air Bleed

33.0025.01



- For Continental whip and Excellence cuspidor
- Use for air bleed
- Snap fit mounting
- 1/8" barb

Brake Actuator Valve

99.0667.00 Actuator valve




- For controlhead brake on A-dec 500 and A-dec 300

Étrangleurs, leviers en plastique

Single Flow Restrictors


025.052.00 3/8" OD tubing



3/8" Plastic

Double Flow Restrictor

23.0173.00



- 1/8" OD tubing
- Used for 3-way air/water syringe flow control
- O-ring added to adjustment screw for shipping. Discard after screw adjustment.

Pressure Indicator

24.0100.03



006.009.00


- Panel mount
- Black and silver bezels included
- Red indicator

Panel Mounting Brackets

026.036.00 Single

42.0701.00 Double

56.0339.00 Double for QDs



Knob, 1/4" ID

027.008.00



Setscrew 007.001.00

- 1/4" ID with set screw
- Replaces all 1/4" knobs
- Aluminum
- Requires 6-32 hex key
- Easy-to-clean design

Toggle, on/off

33.0031.02 pin and toggle kit, Surf 4

22.0462.02 Toggle, Surf 4


22.0462.03 Toggle, White

011.038.01 Pin (5)





Toggle, Momentary

1,5	33.0037.02 Toggle and Pin Kit	
1	011.038.01 Straight pin, .062 Dia x .430, (5)	
2	29.0830.00 Stem w/O-ring, 2 way	
3	006.009.00 Hex nut 15/32-32" x 9/16" x 3/32"	
4	006.069.00 Knurl nut, 15/32-32"	
5	Momentary toggle lever	
	33.0036.00, Black	
	33.0036.01, Gray	
	33.0036.02, Surf	
	33.0036.03, Gray 3	



Canules et embouts seringues

N'utiliser que des embouts A-dec



Standard Syringe Tip

23.0872.01 (5)


23.0872.02 (25)



- Autoclavable
- Standard size 3¹/₄" (83 mm)
- Specially designed for A-dec syringes

Syringe Tip Retainer Non-Locking Kit


23.1112.00 Non-locking kit



- Autoclavable
- U.S. Patent #5,433,485

Syringe Tip Retainer Locking Kit


23.1090.00 Locking kit



- Autoclavable
- Embout à utiliser pour fonction de rétracteur de joues. Il est bloquant.


Autoclavables traditionnelles

Autoclavable Traditional Syringe (beginning October 1995)




Autoclavable Traditional Syringe (prior to October 1995)

Service Parts Only




Terminal Barb Assembly Non-Quick-disconnect

23.1110.00 2-barb



Terminal Barb Assemblies Quick-disconnect

23.1067.00 2-barb




Standard Syringe Kits

	SYRINGE ASSEMBLIES (Include syringe, tubing and flow control clamp)	PART NUMBER SURF	PART NUMBER GRAY 3	PRICE
	With Terminal Barb Assembly 7' (2,134 mm)	23.1150.00	23.1302.00	
	With QD Terminal Barb Assembly 7' (2,134 mm)	23.1195.00	23.1303.00	

Autoclavables continentales


Autoclavable Continental Syringe

(Beginning October 1995)




Terminal Barb Assembly Non-Quick-disconnect

23.1110.00 2-barb



Terminal Barb Assemblies Quick-disconnect

23.1067.00 2-barb



Autoclavable Continental Syringe

	SYRINGE ASSEMBLIES (Include syringe and tubing)	PART NUMBER SURF	PART NUMBER GRAY 3	PRICE
	With Terminal Barb Assembly 84" (2,134 mm)	23.1170.00	23.1304.00	
	With QD Terminal Barb Assembly 84" (2,134 mm)	23.1247.00	23.1305.00	

Boutons et kits de réparation

Button Valves


23.1232.01 Buttons (2)



- Requires spring 013.064.01 (10)

Autoclavable Syringe Valve Kit


23.1234.00



- For syringes beginning October 1995
- Includes button valves, springs, retaining screw, hex key, and silicone lubricant

2-Button Conversion Kit

23.1194.00



- Converts from single autoclavable button to two autoclavable buttons
- Includes: (2) valve assemblies, (1) spacer, screw, (2) springs, silicone lubricant, and hex key

Autoclavable QD Syringe Service Kit

23.1235.01



- Includes syringe tip, screw, pin, syringe nut with spacer, springs, button valves, stem assemblies, terminal barbs, silicone lubricant and O-rings

Kits de réparation

Button Valve

23.0020.02 (2)



Syringe Button Diaphragms


23.0220.02 (10)



- For use on 23.0020.02 syringe button assy

Uni-Flo and Tri-Flo Syringe Valve Core Service Kit


23.0018.00



- Includes four valve cores, two syringe button diaphragms, and valve core removal tool

Uni-Flo and Tri-Flo Syringe Service Kit


23.0019.00



- Includes two syringe tips, nut, four valves, two diaphragms, two O-rings, valve core removal tool, and lubricant

Uni-Flo Syringe Service Kit


23.0027.00



- Includes one button assembly, two valves, two O-rings, valve core removal tool, and lubricant

A-dec syringe tool


98.0010.00



A-dec's syringe tool incorporates:

- A valve core removal tool
- A 3/32" hex wrench for syringe flow adjustment
- A 5/32" hex wrench for removing and installing the syringe tip insert

Composants et tuyaux

Syringe Tubing		PART NUMBER SURF	PART NUMBER GRAY 3	
	Straight 3-hole 7' (2,134 mm)	024.155.02	024.223.02	
	Straight 3-hole 12' (3,657.6 mm)	024.155.12	024.223.12	
	Assembly 7' (2,134 mm) (Includes collet, bushing, retainer, retainer spring, and uni-clamps)	23.1208.00	23.1208.01	

Replacement Handle

23.1015.00



- Use for:
Autoclavable, and Continental
autoclavable syringes

Syringe Strain Relief Repair Kit

90.0563.00 Repair kit

010.073.01 Spring (5)



- Includes (1) O-ring and
(2) retainer springs

Pédales montées




A-dec 500 Wet/Dry Foot Control with Chip Blower/Accessory Button

A-dec 500 wet/dry foot control

with chip blower/accessory button

- Modulates variable speed operation with foot pressure applied in any position
- Use for air turbine, electric handpiece, A-dec camera, as well as ultrasonic instruments
- Dark Surf or Gray 3
- For use with A-dec 500 delivery systems only

PRODUCT DESCRIPTION	PART NUMBER SURF	PART NUMBER GRAY 3	PRICE
A-dec 532/533 foot control 5' (1,524 mm)	38.1805.00	38.1829.00	
Tubing replacement kit	90.1104.00	90.1104.01	
A-dec 541/542 foot control, 16' (4,876 mm)	38.1805.01	38.1829.01	




A-dec 300 Wet/Dry Foot Control

A-dec 300 wet/dry foot control

- Variable speed operation with foot pressure applied in any position
- Toggle operates On/Off control for coolant water
- Use for both the air turbine or electric handpiece

PRODUCT DESCRIPTION	PART NUMBER	PRICE
A-dec 332/333 foot control 5' (1,524mm) Gray 3	38.1834.00	
A-dec 334/335 foot control 9' (2,743mm) Gray 3	38.1821.00	
A-dec 300 euro conversion kit	38.0825.00	




A-dec 300 Wet/Dry Foot Control with Chip Blower

A-dec 300 wet/dry foot control

with Chip Blower

- Variable speed operation with foot pressure applied in any position
- Toggle operates On/Off control for coolant water
- Use for air turbine, electric handpiece, A-dec camera, as well as ultrasonic instruments

PRODUCT DESCRIPTION	PART NUMBER	PRICE
A-dec 332/333 foot control 5' (1,524mm) Gray 3	38.1862.00	
A-dec 334/335 foot control 9' (2,743mm) Gray 3	38.1838.00	
A-dec 300 euro conversion kit	38.0825.00	




Foot Control

Foot control III

- Variable speed operation with foot pressure applied in any position
- For use with Cascade or Radius units equipped with the Master Series touchpad option
- Standard tubing: 7' (2,134 mm) Dark Surf

PRODUCT DESCRIPTION	PART NUMBER	PRICE
Cascade foot control	38.1700.00	




Wet/Dry Foot Control

Wet/dry foot control III

- Variable speed operation with foot pressure applied in any position
- Toggle operates On/Off control for coolant water
- Standard tubing: 7' (2,134 mm) Dark Surf
- Use for both the air turbine or electric handpiece

PRODUCT DESCRIPTION	PART NUMBER	PRICE
Dark Surf	38.1603.00	
Euro wet/rdy, Dark Surf	38.1603.01	
Foot control III euro conversion kit	38.1764.00	



Wet/Dry Foot Control with Chip Blower/Accessory Button

Wet/dry foot control III

with chip blower/accessory button


- Variable speed operation with foot pressure applied in any position
- Toggle operates On/Off control for coolant water
- Button operates chip blower valve or air-actuated electric switch for ultrasonic instruments
- Standard tubing: 7' (2,134 mm) Dark Surf

PRODUCT DESCRIPTION	PART NUMBER	PRICE
Chip blower, Dark Surf	38.1601.00	
Euro wet/dry/chip blower, Dark Surf	38.1601.01	
Foot control III euro conversion kit	38.1764.00	

Composants des pédales et kits de reconstruction

Foot Control I Valve Service Kit


90.0010.00



- For use with Foot Control I Valve (38.0659.00)

Foot Control III Valve


38.0761.00



- With relay
- Replaces 38.0607.01 and 38.0659.00
- Not for A-dec 500 foot control

Lever Foot Control Service Kit


90.0313.00



- For use with Lever Foot Control Valve (38.0201.00)

Foot Control II Valve Service Kit


90.0312.00



- For use with Foot Control II Valve (38.0607.01)

Foot Control III Valve Service Kit

90.0593.00



- For use with Foot Control III Valve (38.0761.00)
- Also use with A-dec 500, 400, and 300 foot controls

Poppet
22.0060.00

Service Tip

FUITE D'EAU SUR LA PIÈCE À MAIN APRÈS AVOIR RELÂCHÉ LA PÉDALE :


Certaines fuites d'eau peuvent être causées par un tuyau de pédale pincé qui ralentit l'échappement de l'air qui bloque le signal air de commande de l'eau.

Composants des pédales et kits de reconstruction

Wet/Dry Toggle (3-Way) Valve


38.0604.00 Surf 6

38.1832.00 Gray 3



Chip Blower Button with 2-Way Valve

38.0610.00 Surf 6




- For use with Foot Controls I, II, and III as a chip blower button

Chip Blower/Scaler/Accessory Button with 3-Way Valve

38.0612.00 Surf 6

38.0612.01 Gray 3




- Use with A-dec 500 foot controls as chip blower/Accessory button
- Use with Foot Controls I, II and III as a chip blower button

Micro Valve

33.0134.00 2-Way


33.0138.00 3-Way



- 2-way (brass ball) for chip air
- 3-way (stainless steel ball) for wet/dry or scaler

Foot Control I and II One-Hole (right) Housing


38.0320.02 Surf 6



- Use with chip blower, ultrasonic instrument, wet/dry or foot-activated water dispensing system

Foot Control III One-Hole (left) Housing

38.0763.00 Surf 6




- Use with Foot Control III with wet/dry toggle

Foot Control Housing, Two-Hole

38.1802.00 A-dec 500, Surf 6

38.0321.02 Surf 6 (FC I, II, III)

38.1802.01 A-dec 500, Gray 3



- Use with wet/dry with chip blower, wet/dry with ultrasonic instrument, or ultrasonic instrument with chip blower


Composants des pédales et kits de reconstruction

Toggle Kit

38.0075.03 Surf 6

38.0075.04 Gray 3

38.0075.00 Black



- Includes pin, toggle, and retaining spring

Rubber Bumper Pad (Feet)


017.007.01 (25)



- 1/2" diameter bumper with self-adhesive back

Foot Control II and III Retaining Rings


38.0237.00 Black



- For Foot controls II and III
- Includes 4 screws 6 x 1/2"

Foot Control Cover


22.0110.00



- Fits all A-dec foot controls
- Chrome plated

Foot Control I Retaining Ring

22.0120.00 Gray




- Includes 4 screws 6 x 1/2"

Valve Assembly Screws

003.078.00 Socket Head Screw (1)

4-40 x 1/4"

- Attaches valve assembly to foot control base

Foot Control Tubing	LENGTH	NO SHEATH	DARK SURF	GRAY 3	
 Tubing Without Sheath	10' (3m)	024.151.01	-	-	
	50' (15m)	024.151.03	-	-	
 Foot Control III Tubing	10' (3m)	-	024.179.01	-	
	50' (15m)	-	024.179.03	-	
 A-dec 500 Foot Control Tubing	10' (3m)	-	024.202.01	024.229.01	
	50' (15m)	-	024.202.02	024.229.02	
 A-dec 300 Foot Control Tubing	10' (3m)	-	-	024.214.01	
	50' (15m)	-	-	024.214.02	

A-dec 500 et 300 cartouches, clapets, diaphragmes

Water Cartridge

38.1780.00 Water cartridge



- For converting from a dry control block position to a wet block
- Includes washer
- For A-dec 500 and 300 control blocks only

Dry Cartridge

38.1783.00 Dry cartridge



- For converting from a wet control block position to a dry block
- Includes washer
- For A-dec 500 and 300 control blocks only

Control Adjustment Stems For A-dec 500

38.1785.01 Drive air stem

38.0526.00 Air coolant stem

38.0516.00 Water coolant stem



Drive air stem



Air coolant stem



Water coolant stem

Flow Adjustment Key For A-dec 500

38.1808.01 Flow adj. key (4)



- Key with O-ring adjusts air and water coolant

Control Block Diaphragm

38.1787.01 Diaphragm, 4-position blk (5)



- For A-dec 500 and 300 control blocks

Control Adjustment Stems for A-dec 300

38.1785.01 Drive air stem

38.1825.00 Air coolant stem

38.1824.00 Water coolant stem



Drive air stem



Air coolant stem



Water coolant stem

Control Block QD Plug

026.160.01 QD, plug (5)



- For A-dec 500 and 300 control block
- Includes O-ring

Accessory Knob Kit

90.1089.00 Accessory knob kit



- Use for accessory item control

Flow Adjustment Key For A-dec 300

38.1823.00




- 7/16" dia. white 2

Bloc Century plus

Century Plus Control Block

38.0540.00




- With pigtails

Century Plus Cartridges

38.0517.00 Air bleed

38.0518.00 Check valve

38.0520.00 Water valve



Air bleed
(brass base)



Check valve
(blue base
with duckbill)



Water valve
(red base)

Drive Air Flow Adjustment Stem

38.0510.00 Stem

035.034.01 O-ring (10)




- Includes O-ring

Water Coolant Flow Adjustment Stem

38.0516.00 Stem


035.034.01 O-ring (10)



- Includes O-ring

Syringe Flow Control Barb

38.0555.00



- For air and water
- 1/8" barb
- 10-32 threaded fitting


Flush Valve Stem

38.0525.00 Stem

035.033.01 O-ring (10)

034.001.01 O-ring (10)

013.072.00 Spring



- Includes O-rings

Air Coolant Flow Adjustment Stem

38.0526.00 Stem

030.003.02 O-ring (10)



- Includes O-ring

Handpiece Control Adjustment Keys

98.0867.02 (5)




- Autoclavable adjustment keys
- For use in Cascade and Decade 4400 handpiece controls

Bloc Century plus

Traditional-Style Century Plus Handpiece Control Kits

38.0631.00 Unitized

38.0375.00 Individual w/holder




- For Cascade handpiece position
- Includes control block, valve, tie bolts, gasket, and hardware

Continental-Style Century Plus Handpiece Control Kits

38.0638.00 Continental

39.1054.00 Whip only




- For Cascade handpiece position
- Includes control block, valve, tie bolts, whip assembly, gasket, and hardware

Pins For Continental Whip

39.1050.00 Short pin

39.1051.00 Long pin




- Short pin used on outer sides of end whips
- Long pin used between whips

Tie Bolt Kit

38.0504.07 3 blocks

38.0504.08 4 blocks


38.0504.09 5 blocks



- Includes (2) bolts, (2) washers, and (2) nuts

Dry Block Conversion Kit


38.0648.00



- Black base
- Contains (2) plugs (one water, one check valve)

Century Plus Control Block Service Kit

38.0537.01 Service Kit



- Includes cartridges, gaskets, stem, diaphragms, cap assembly, actuator, washers, spring, and screws
- Order IOLS diaphragm for systems that also include a Fiber-optic Light Source (75.0911.01)

Continental Whip Spring Washer Kit

90.0446.00 (10)



- Concave side of washer faces the whip

Continental Whip Tension Springs

013.015.00 Red 3 lb. spring


013.027.00 Green 4 lb. spring



- Red 3 lb. spring is standard on Cascade
- Green 4 lb. spring used for greater return pull on whip arm

Radius Delivery System Coolant Extension Knob


83.0217.00 Coolant Extension Knob



- Use if existing knobs are too short

Century Plus Manifold Assembly for Cascade


38.0524.00 for Cascade



Bloc Century II

Century II Control Block


38.0223.02 Block assembly
38.0234.02 Kit



- Block assembly includes control block, (2) knobs (one gray, one black), and gasket
- Kit includes control block, valve, tie bolts, gasket, holder, and hardware

Water Valve Cover and Stem Assembly


29.0831.00 Stem w/O-rings
38.0181.00 Valve cover
38.0054.02 Diaphragm (10)



- For Century II control blocks

Century II Control Block Repair Kit


90.0321.00 Kit



- Includes diaphragms, gaskets, O-rings, stem spring, and screws

Push On Knob Kits for Coolant Stems

027.017.02 Blue dot for water
027.018.02 Yellow dot for air



- Each kit includes both gray and black knobs w/color dot
- Choose yellow for air coolant and blue for water coolant

IC Control Block


38.0467.00 Block assembly



- Block assembly includes control block, barbs, and caps

Valve Cover Kit


90.0304.00 Kit
38.0178.00 Cover w/screws
38.0179.02 Diaphragm (10)



- Includes valve cover, diaphragm, O-rings, and screws

Tie Bolt Kit


38.0193.01 3 blocks
38.0194.01 4 blocks



- Includes (2) bolts, (2) washers, and (2) nuts

Air and Water Coolant Stem


33.0103.00 Stem
030.004.02 O-ring (10)



- Includes O-ring (030.004.02)

Compression (Stem) Spring


013.025.00



Bloc Century II

3-Way Cartridge


33.0100.02 Surf 4 toggle



- Use for Master On/Off and coolant water air signal On/Off functions in Century II master block

Cartridge Needle Valve Kit


33.0101.02



- Includes (2) knobs (one gray, one black)
- Use for coolant air metering in Century II master block


Adjustable Screw w/O-ring

29.0100.00



Air Bleed Restrictor


38.0316.00 Restrictor
030.004.02 O-ring (10)



- Includes O-ring

Flush/Non-Retractable Handpiece Conversion Kits


14.0326.00 Non-retraction kit



- Add flush capability to handpiece tubing
- 1 handpiece kit includes teed check valve, uniclamps, and instructions

On/Off Indicator


24.0105.01 Silver bezel



- Used for On/Off indicator in Century II master block

Key Adapter with Screw


39.0820.00



- For use on IC control heads to adjust scaler blocks

Air Bleed Filter


19.0001.03 (2)



- 10 microns

Air Bleed Restrictor Barbs

023.036.00 .005
023.070.00 .013
004.005.02 Washer (10)




- Requires washer (004.005.01)

Bloc Performer

Performer Control Block

38.0729.00 Block w/pigtails


38.1775.00 Service kit, Perf



- Service kit includes: diaphragm, water relay assembly, and water valve actuator

Performer Control Block Diaphragm

38.0711.01 (5)



Water Relay Valve Cartridge and Actuator

38.0717.00 Cartridge

38.0718.00 Actuator




Cartridge



Actuator

Performer Knob, Push On

027.057.00




Coolant Flow Adjustment Stems


38.0712.00 Water stem

38.0713.00 Air stem

030.004.02 O-ring (10)




- Water stem
Includes O-ring



- Air stem
Includes O-ring

Performer Pressure Relief Kit

90.1020.00 Kit



- Includes muffler assembly, valve assembly, (5) uni-clamps, and installation instructions

Flow Control Screw

38.0766.02 Assembly (5)

030.004.02 O-ring (10)



- Assembly includes screw with O-ring

Service Tip

Use a Performer pressure relief kit (90.1020.00) to prevent trapped air and water from exhausting through the high speed handpiece when it is lifted from the holder.

Fouets, poignées et supports

12 O'Clock Doctor's Holder with Valve

- 99.0711.00 541 Dr.'s holder, Surf 4
- 99.0711.01 541 Dr.'s holder, White 2



- For Doctor's handpieces on a 12 o'clock delivery system

Whip Assembly Kit for A-dec 500

- 77.0291.01 Whip kit, Surf 4
- 77.0291.02 Whip kit, White 2
- 013.123.00 Spring, whip
- 77.0063.00 Whip cover, Surf 4
- 77.0063.01 Whip cover, White 2



- Includes whip assembly, cover, and spring

Whip Components for A-dec 300

- 77.0769.00 Whip hook
- 90.1205.00 Whip kit A-dec 300



- Includes whip assembly, cover, and spring


A-dec 500 Accessory Holders

- 99.0681.00 Holder, Surf 4
- 99.0681.01 Holder, White 2
- 99.0684.00 Holder, auto-electric, Surf 4
- 99.0684.01 Holder, auto-electric, White 2
- 99.0708.01 12 O'clock Holder, White 2




Pompe à salive

- 99.0653.01 Holder, Surf 4
- 99.0653.03 Holder, White 2
- 99.0673.01 Kit, holder, auto-electric, Surf 4
- 99.0739.01 Kit, holder, auto-electric, White 2



Grosse canule 15 mm HVE Holder

- 99.0675.01 Holder, auto-electric, Surf 4
- 99.0738.01 Holder, auto-electric, White 2



Poignées 300 & 500

- 77.0104.03 500 Brake Handle Assembly, Grey
- 77.0104.01 500 Brake Handle Assembly, Surf
- 90.1281.00 300 Brake Handle Assembly



Flex Holder Assembly

- 99.0705.00 Syringe/SE holder, Surf 4
- 99.0705.01 Syringe/SE holder, White 2
- 99.0707.01 15 mm HVE holder, White 2




Supports pièces à main et collecteur

3-Position Cascade Unitized Holders

99.0619.00 Left

99.0618.00 Right




- Includes unitized holder (Surf 4), air bleed valves, tubing, and necessary hardware

Plugs for Unitized Holder


99.0607.00 Hole plug, Surf 4

99.0587.00 Slot plug, Surf 4

75.0083.00 Continental plug



Hole Plug




Slot Plug



Continental Plug

Cascade Individual Holder

99.0583.00




45.0403.00
Friction pad

- Fits 1/2" bar
- Includes individual Surf 4 holder, air bleed valve, tubing friction pad, set-screw, and necessary hardware

Auto Holders

99.0537.01 Gray assembly

99.0541.01 Gray holder only




- Fits 1/2" bar
- Includes tubing and necessary hardware

Handpiece Oil Collection System

24.0410.00

90.1110.00 A-dec 500 Surf


90.1110.01 A-dec 500 & A-dec 300 White



- Collects oil expelled from the handpiece exhaust
- Accommodates up to three handpiece tubing assemblies

Individual Cascade Assistant's Holder

99.0584.00 For syringe and SE




45.0403.00
Friction pad

- Fits 1/2" bar
- Includes individual Surf 4 holder, friction pad, setscrew, and necessary hardware

Valves de pré-sélection

2-way Air Bleed Valves Use for Air bleed		PART NUMBER	
	for Cascade Unitized Holder 1/8" barb, snap fit mounting	33.0132.00	
	for Cascade Individual Holder and Excellence Holder 1/8" barb, snap fit mounting	33.0025.00	
	for Cascade Continental Whip Holder 1/8" barb, snap fit mounting	33.0025.01	

Actuators		PART NUMBER	
	for Cascade Unitized Holder	99.0590.00	
	Actuator with Valve Assembly for Cascade Unitized Holder	99.0592.00	

2- or 3-way Micro Valves		PART NUMBER	
	2-way, brass ball 1/8" barb	33.0134.00	
	3-way, stainless steel ball 1/8" barb	33.0138.00	

Holder Valves		PART NUMBER	
• replace valve by detaching tubing from control block			
	Performer	99.0724.00	
	A-dec 500, Surf	99.0654.02	
	A-dec 500 & A-dec 300, White	99.0654.03	
	A-dec 541, Surf	99.0696.01	
	A-dec 541, White	99.1135.01	

A-dec 500 et 300

PCB Standoffs

90.1082.00 Standoff (5)



PCB Jumpers


041.427.01 Jumpers (10)



Deluxe Touchpad Dual Solenoid Manifold Assembly for A-dec 500

43.0057.00 Manifold, dual

041.659.01 Solenoid, only



- Assembly includes two 5 VDC solenoid valves and dual manifold

A-dec 571 Dental Light Intensity Switch

43.0054.00 Switch assembly



Data Cable Adapter

041.684.00 Adapter



- Data cable adapter
- 6-pin conductor, one female to two females

Wago Connector

90.1176.00 Wire connectors (6)



- For in-line connections

A-dec 500 et 300 câbles RJ, électrovanne et module relais

A-dec 500 and 300 Data Cables and Data Board	PART NUMBER	PRICE
6 Pos, 1.5' (457 mm)	041.637.00	
6 Pos, 3.5' (1,029 mm)	041.637.06	
6 Pos, 7' (2,134 mm)	041.637.04	
6 Pos, 8.5' (2,589 mm)	041.637.05	
6 Pos, 13' (4,010 mm)	041.637.02	
6 Pos, 25' (7,614 mm)	041.555.00	
Kit, 25' data cable and 2 adapters	041.554.00	

Printed Circuit Board Kits	PART NUMBER	PRICE
A-dec 500		
Chair	90.1072.00	
Control head	90.1076.00	
Cuspidor	90.1079.00	
Dental light	28.1577.00	
A-dec 300		
Chair	90.1177.00	
Quad volt intraoral light source	90.1168.00	
Cuspidor	90.1079.00	

A-dec 311 & 411 Chair Solenoid Assembly

90.1225.00 Solenoid kit



- 21.6 VDC

A-dec Relay Module


90.1080.00 Kit, A-dec relay



UTILISER CE MODULE SI :

- commande de moteur d'aspiration
- commande d'éclairage ambiant
- commande caméra intra-buccale

Platines électriques, 25 W, 80 W et 300 W



300-watt power supply

- Required for intraoral camera and electric motors
- Provides power for multiple electrical accessories such as fiber-optics, low voltage water heater/heated syringe tubing, dental light, ultrasonic instrument, and curing light.
- Size – Fits the large, contoured, and recessed floor boxes
- QD Connectors
- Maintenance – Uses resettable circuit breakers rather than fuses

PRODUCT DESCRIPTION	PART NUMBER	PRICE
25 watt, 220-240V, kit	28.1481.00	
80 watt, 220-240V	47.2031.00	

PRODUCT DESCRIPTION	PART NUMBER	PRICE
A-dec 500 Chair Power Supplies		
511 300 watt, 220-240V	90.1075.00	
Cascade/Decade Chair Power Supplies		
Cascade 300 watt, 240V	28.1436.00	
Cable adapter kit (used to retrofit from 300 watt power supply manufactured before June 1998)	90.1012.00	
NOTE: When purchasing a power supply, consider both present and future electrical accessories. Buying a 300 watt power supply now will make it economical to add electrical accessories in the future.		

Câbles d'alimentation pour 541, 542 et 545

Remote Power Supply cables

Power supply cables are required if you are connecting a delivery system (A-dec 541, 542) or assistant's instrumentation (A-dec 545) to a remote power supply



ORDERING CABLES

- Refer to the table below to determine your cable needs
- See the Ordering Information Table for cable descriptions and prices

POWER SUPPLY LOCATION	DELIVERY SYSTEM/ASSISTANT'S INSTRUMENTATION			
	541	542	545	545 WITH 542
300-watt, 511 Chair	77.0707.00	77.1239.00	77.0708.00	(1) 77.0708.00 (1) 77.1239.00
300-watt, 311 Chair	77.0980.00	77.1240.00	77.0981.00	(1) 77.0981.00 (1) 77.1240.00
300-watt, remote	(1) 43.0328.01 (1) 43.0093.00	43.0329.00	43.0095.00	(1) 43.0329.00 (1) 43.0093.00
1/8" Air Line Extension Kit (25' [7,620 mm])	90.1129.00			

ORDERING INFORMATION / PRODUCT DESCRIPTION	PART NUMBER	PRICE
Air Line Kit, remote power supply	90.1129.00	
Cable, 15' (4,572 mm); 6-pin to 6-pin; Assistant's	43.0093.00	
Cable, 15' (4,572 mm); 6-pin to 4-pin; Assistant's	43.0095.00	
Cable, 15' (4,572 mm); 4-pin to 4-pin; Delivery	43.0328.01	
Cable, 15' (4,572 mm); 6-pin to 6-pin; Delivery	43.0329.00	
545 to 511 Power Cable Kit	77.0708.00	
545 to 311 Power Cable Kit	77.0981.00	
542 to 511 Power Cable Kit	77.1239.00	
542 to 311 Power Cable Kit	77.1240.00	
541 to 511 Power Cable Kit	77.0707.00	
541 to 311 Power Cable Kit	77.0980.00	

Quad Voltage Intraoral Light Source



IOLS kit for A-dec 300

- Fiber-optic intraoral light source
- Requires fiber-optic handpiece tubing, 24vac
- All IOLS come equipped to light (2) handpieces
- A-dec 300 kit equipped to light (4) handpieces
- To add additional handpieces order the QVIOLS Switch Kit

PRODUCT DESCRIPTION	PART NUMBER	PRICE
Quad Volt IOLS kit for A-dec 300 (equipped to light (4) handpieces)	90.1168.00	

Claviers



A-dec 300/500
Standard Touchpad



A-dec 300
Deluxe Touchpad



A-dec 500
Deluxe Touchpad

A-dec 300 and 500 Dessus clavier seul	PART NUMBER SURF	PART NUMBER WHITE	PRICE
Standard Touchpad Cover Kit	90.1439.00	90.1440.00	
A-dec 500 Deluxe Touchpad Cover Kit	90.1441.00	90.1442.00	
A-dec 300 Deluxe Touchpad cover Kit	–	90.1443.00	
Disposable Wrist Strap	041.806.00		

A-dec 300 Touchpads

- For use with A-dec 300 delivery systems only

PRODUCT DESCRIPTION	PART NUMBER	PRICE
Complet		
Standard touchpad kit, 334/335	90.1196.00	
Deluxe touchpad kit, 334/335	90.1197.00	
Standard Touchpad Kit, 352/353	90.1184.00	
Dessus clavier et carte électronique seulement		
Standard touchpad kit, 334/335, 352/353	90.1078.01	
Deluxe touchpad kit, 334/335	90.1185.00	

A-dec 500 Touchpad Upgrade Kits

- For use with A-dec 500 delivery systems only
- Replace an existing touchpad
- Add a new touchpad


PRODUCT DESCRIPTION	PART NUMBER SURF	PART NUMBER WHITE	PRICE
Complet			
Standard Touchpad, 542	77.0297.00	77.0297.01	
Standard Doctor's Touchpad, 541	90.1127.00	90.1127.01	
Standard Assistant's Touchpad, 541/545	90.1128.00	90.1128.01	
Deluxe Touchpad Upgrade Kit, 532/533/542	90.1133.00	90.1133.01	
Deluxe Touchpad Upgrade Kit, 541	90.1134.00	90.1134.01	
Dessus clavier et carte électronique seulement			
Standard Touchpad	90.1078.00	90.1078.01	
Deluxe Touchpad	90.1135.00	90.1135.01	

Interrupteurs électriques

Air Electric Switch

(Normally Closed Or Open)

044.170.00




- 40 psi
- For 25 W, 150 W, and 300 W, power supplies
- For master control, dental light
- Rating: 15A, 100-240VAC
- For A-dec 300/400/500 power supplies

Air Electric Switch

(Normally Open)

40.0250.00 Kit




- 10 psi
- For scaler
- Rating: 5A, 24V AC/DC
- Kit includes switch, 1/8" barbed tee, 5' clear tubing uni-clamps, U-clip wire nuts, terminal QDs check valve

Air Electric Switch

(Normally Closed)

044.159.00 Switch only




- 10 psi
- For Focus fiber-optic and Dual Voltage intraoral Light Source
- Rating: 5A resistive, 24V AC/DC
- Kit includes: switch, 1/8" barbed tee, 5' clear tubing, uni-clamps, U-clip, wire nuts, terminal QDs

Pompes à salive et grosse canule, adaptateur


Grosse canule 11 mm

11.1335.00 without tailpiece




Pompe à salive

12.1258.00 without tailpiece




Grosse canule 15 mm

12.1257.00 without tailpiece



Tubing Adapter $\frac{5}{8}$ " to $\frac{1}{2}$ " pour 11 mm


40.0145.00



Écouvillons, entonnoir, canules chir. et kits de conversion

Pour simulateurs


- 10.0014.01 Tip (5)
- 10.0015.00 Tip assembly



- Stainless steel
- 7/16" diameter OD

Surgical Aspirator Tips


- 11.0070.00 1/8" OD
- 11.0090.00 3/16" OD
- 11.0170.00 5/16" OD (Ohio)



- Autoclavable stainless steel tip
- Adapts to the AVS and central vacuum HVE valve, includes cleaning brush


Surgical Tip Brushes

- 049.005.01 1/8" (6)
- 049.002.01 3/16" (6)



Vacuum Valve Brushes


- 049.001.01 1/2" (6)
- 049.004.01 3/4" (6)



- Use 1/2" for AVS
- Use 3/4" for 15 mm European and standard

Oral Cups and Liners pour 11 mm

- 11.0440.00 Molded
- 11.0490.00 Liners (1000)




- Molded and metal oral cups require paper liners

A-dec 500




DESCRIPTION	SURF	WHITE	PRICE
551 Vacuum Tower Assembly	90.1148.00	90.1148.01	


DESCRIPTION	SURF	WHITE	PRICE
12 O'Clock Vacuum Canister Assembly	90.1284.00	90.1284.01	



REF	DESCRIPTION	SURF	WHITE	PRICE
1	Vacuum Tower Cap includes (1) O-ring	77.0254.01	77.0254.03	
2	O-ring (10)	030.035.02		
3	Screen, trap (10)	11.1284.01		
4	Triple Vacuum Connector includes (3) O-rings	77.0453.00	77.0453.01	
5	O-ring Connector Kit includes (3) O-rings	90.1357.00		
6	15 mm Vacuum Cap includes (1) O-ring	11.1280.01	11.1280.02	
7	551 Vacuum Tower Base w/Funnel includes (1) red O-ring	77.0362.01	77.0362.03	
8	O-ring, Red (10)	035.053.01		
9	12 O'Clock Vacuum Canister Base	77.1104.00	77.1104.01	
10	551 Vacuum Pipe includes (1) red O-ring	77.0121.01	77.0121.03	




A-dec 300




DESCRIPTION	SURF	WHITE	PRICE
300 Vacuum Canister Assembly	-	77.0941.01	



DESCRIPTION	SURF	WHITE	PRICE
351 Vacuum Canister Assembly	-	90.1280.00	




REF	DESCRIPTION	SURF	WHITE	PRICE
1	Vacuum Tower Cap, includes (1) O-ring	77.0254.01	77.0254.03	
2	O-ring (10)	030.035.02		
3	Screen, trap (10)	11.1284.01		
4	Triple Vacuum Connector, includes (3) O-rings	77.0453.00	77.0453.01	
5	O-ring Connector Kit, includes (3) O-rings	90.1357.00		
6	15 mm Vacuum Cap, includes (1) O-ring	11.1280.01	11.1280.02	
7	300 Vacuum Canister Base		77.0748.00	
8	300 Vacuum Canister Blank Cover		77.0724.00	
9	351 Vacuum Canister Base, includes (1) O-ring		77.1047.00	



Collecteurs


Cascade Solids Collector Single HVE

For Std or 15 mm HVE




Pinnacle Screen Kit

11.1194.00 For single HVE



- Includes (3) screens and replacement cap 11.1192.00 for single HVE
- Pinnacle screens are available from Pinnacle; order #5500-1


Cascade Solids Collector 15 mm HVE with Pinnacle Screen Kit



Central Vacuum Solids Collector

17.0600.02 Canister assy

028.008.00 Plug/Cap, 1/2" x 1/4"



A-dec 500, 400 et 300

Solenoid, Cuspidor

041.660.00 Solenoid w/O-ring
030.016.02 O-ring (10)



A-dec 500 and 400 Bowl Rinse and Cup Fill Spouts

77.0042.01 Spout, bowl rinse, Surf
77.0042.04 Spout, bowl rinse, Gray 3
77.0043.01 Spout, cup fill, Surf
77.0043.04 Spout, cup fill, Gray 3



A-dec 500 Cuspidor and Post Trim Rings

77.0098.00 Cuspidor, Surf 4
77.0098.01 Cuspidor, White 2
77.0095.00 Post, Surf 4
77.0095.01 Post, White 2



Cuspidor



Post

Bowl Rinse/Cup Fill/Safety Stop Switch Assembly

43.0010.00 Switch assembly



• Operates the cuspidor bowl rinse and cup fill options


A-dec 300/400/500 Spout Gaskets

77.0044.00 Surf 4
77.0044.01 Gray 3



A-dec 300/400/500 Cuspidor Printed Circuit Board Kit

90.1079.00



A-dec 300 Cuspidor Bowl

77.0730.01 Bowl with O-rings
030.224.01 O-ring (10)



A-dec 300/400/500 Cuspidor Bowl Screens

75.0035.03 Bowl screen, White 2 (5)
75.0035.01 Bowl screen, Surf (5)



A-dec 300 Bowl rinse and Cup Fill Spouts and cup landing

77.1033.00 Spout, bowl rinse, with O-ring
77.1034.00 Spout, cupfill, with O-ring
90.1213.00 Cup landing




P-trap for a-dec 300 gravity drain

77.0923.00




Cascade

Cascade Vitreous China Cuspidor Bowl



REF NUMBER	DESCRIPTION	PART NUMBER	PRICE
1	Contoured cup filler water spout	75.0039.00	
2	Spout gasket (5)	75.0145.01	
3	O-ring, .301 ID x .070 W (10)	030.011.02	
4	Valve assembly, cuspidor bowl rinse	12.1028.00	
5	Bowl rinse spout	12.1031.00	
6	Bowl screen (standard) (5)	75.0035.01	
7	Expansion mounting nut (4)	006.133.01	
8	2-way micro valve, bowl rinse	33.0134.00	
9	3-way micro valve, cup fill	33.0138.00	
10	Actuator cover	75.0091.00	


Cascade Vitreous China Cuspidor Base



REF NUMBER	DESCRIPTION	PART NUMBER	PRICE
11	Timed bowl rinse relay valve	12.0913.00	
12	Switch, cuspidor stop	12.0915.00	
13	3-way restricted diaphragm valve (replaces 33.0139.00)	33.0160.00	
14	Cuspidor cup timer stabilizer kit (replaces 12.0954.01 valve)	90.0456.00	
15	Pilot air mini valve	12.0953.00	
16	Drain seal	12.1024.00	
17	Mounting screw assembly (4)	12.1054.02	
18	Drain tube	12.1034.00	
19	Cup fill relay	12.0934.00	


Cascade

Cascade Phenolic Cuspidor Bowl



REF NUMBER	DESCRIPTION	PART NUMBER	PRICE
1	Contoured cup filler water spout	75.0039.00	
2	Bowl screen (standard) (5)	75.0035.01	
3	2-way micro valve, bowl rinse	33.0134.00	
4	3-way micro valve, cup fill	33.0138.00	


Cascade Phenolic Cuspidor Base



REF NUMBER	DESCRIPTION	PART NUMBER	PRICE
7	Timed bowl rinse relay valve	12.0913.00	
8	O-ring, .301 ID x .070 W (10)	030.011.02	
9	Switch, cuspidor stop	12.0915.00	
10	3-way restricted diaphragm valve	33.0160.00	
11	Cuspidor cup timer stabilizer kit	90.0456.00	
12	Pilot air mini valve	12.0953.00	
13	Drain seal	12.1024.00	
14	Cup fill relay	12.0934.00	

Performer

Performer Cuspidor Bowl




REF NUMBER	DESCRIPTION	PART NUMBER	PRICE
1	Cup fill spout	12.0985.00	
2	Bowl rinse spout	12.0986.00	
3	Cuspidor bowl assembly	76.2011.00	
4	Bowl screen (5)	75.0035.01	
5	750 ml water bottle, cap, and supply tube	14.0416.00	
6	Drain seal	12.0991.00	
7	Drain tube	12.0989.00	
8	Hole plug, 5/8"	018.086.00	
9	Fill/rinse tilt rod (2)	12.0977.01	
10	Water actuator housing	12.0987.00	
11	Water spout seal	12.0988.00	
12	Spring, .420 OD x .88	013.004.00	
13	Actuator	12.1013.00	
14	Retaining ring (10)	010.045.02	
15	Diaphragm retainer	12.0983.00	
16	Water manifold diaphragm (4)	12.0982.01	
17	Performer lower cuspidor housing	76.2020.02	

Composants

Bowl Screen

75.0035.01 (5) Surf

75.0035.03 (5) White




- For Cascade, Radius, Performer, A-dec 500, A-dec 400, A-dec 300, and A-dec 200 cuspidors.

Bowl Screen Assembly

42.0085.00 Assembly


42.0090.00 Umbrella only



- Includes replacement umbrella and strainer
- Fits Decade 7283, Decade, Component 7280, Century II, and Excellence cuspidors

8" Cuspidor Bowl


40.0849.00



- 8" glass cuspidor bowl
- Fits Decade 7283, Decade, Component 7280, Century II, and Excellence cuspidors
- Includes (2) O-ring packages (030.228.02)

Cup Holder Assembly


12.0677.01 Surf



- Fits Excellence cuspidor

Cup Filler Tube Assembly


12.0783.00



- Fits Decade 7283 Cuspidor
- Includes (2) O-ring packages (030.011.02)

Bowl Rinse Tube Assembly


12.0784.00



- Fits Decade 7283 Cuspidor
- Includes (2) O-ring packages (030.011.02)

Bowl Rinse Spout

42.0008.00




- Spout with O-rings
- Fits Decade, Component 7280, and Century II post mounted cuspidors
- Includes (2) O-ring packages (030.011.02)

Composants

2-Way Push Button Hex Valve


33.0108.02



- For Decade 7283 Cuspidor cup fill and bowl rinse

Toggle Valve

42.0683.04



- For Decade and Century II
- Surf 4 plastic toggle

006.069.00
\$0.33

006.009.00
\$0.22

Long Stem Air Bleed Valve Assembly

33.0025.01




- For Excellence cuspidor cup fill
- For continental bracket and whip assembly

Actuator Rod Assembly, Cuspidor Tower

12.0709.00 Assembly


030.002.02 O-ring (10)



- Fits Excellence cuspidor

2-Way Micro Valve


33.0134.00



- For Cascade bowl rinse and Excellence
- Brass ball

3-Way Micro Valve


33.0138.00



- For Cascade cup fill and Excellence
- Stainless steel ball

2-Way Valve


38.0086.00



- For vacuum cuspidor
- Balanced

Limit Switch


044.184.01



- For Decade 7283 Cuspidor

Cascade Cuspidor Arm Pivot Hub

12.0931.00



Composants de la LED 570

BARRIÈRES PROTECTRICES

A-dec recommande l'utilisation de barrières protectrices référence 54.0477.01 (pkg. de 10) afin de mieux préserver votre investissement

LED Light Upgrade

28.1925.00 - Surf

28.1924.00 - Grey




- Upgrade for existing A-dec halogen lights
- Compatible with most A-dec lights after 1993
- See note below regarding protective barriers

LED Mode Button


28.1862.00 - Cover only

90.1348.00 - Button kit




LED Indicator Circuit Board Plus Cable

90.1345.00




Light Shield And Gasket Kit

90.1343.00



LED Front Housing

90.1390.00



LED Pivot Shaft Assembly

90.1341.00



Light Handle Kit

90.1342.00 - Surf

90.1342.01 - Grey



LED Neck Pivot Assembly

90.1349.00



CABLE ASSEMBLIES	PART NUMBER	PRICE
LED DCS Dental Light Cable 56"	90.1322.00	
LED DCS Dental Light Cable 82"	90.1322.01	
LED DCS Dental Light Cable 88"	90.1322.02	
LED DCS Dental Light Cable 136"	90.1322.03	
LED DCS Dental Light Cable 220"	90.1322.04	

Composants de l'éclairage 300

Light Bulb

041.709.00 Halogen bulb 95W, 17V



- For A-dec 571, 371, 6300, and pre-Cascade dental lights
- Halogen bulb 95 W, 17V

Flexarm End Cap

75.2018.00 Light bulb holder



- Holder for additional bulb
- Uses 041.709.00 bulb
- For A-dec 371 lights prior to serial number F052936

Bulb Cap

28.1635.00 Bulb cap



Light Shield

28.1634.00 Shield



A-dec 371 Light Vertical adjustment Kit

90.1220.00



A-dec 300 Light post cover Plugs

75.0150.00 hole plug

75.2251.00 Slotted plug for monitor cables




Light Post Mount Cover Plug

- For A-dec 300 intermediate light post

A-dec 300 Intermediate Light Post

28.1624.00 Intermediate post



Trim Ring, A-dec 300 Light Post


77.0876.00 Trim ring



A-dec 500, cascade, 6300 et Performer

Cascade and 571 Dental Light Bulb


- 75.0084.00 Holder Surf 4
- 75.0084.01 Holder White 2
- 28.1004.00 Bulb and holder Surf 4
- 28.1004.01 Bulb and holder White 2



- For Cascade and A-dec 571 dental lights
- Includes holder and halogen bulb 95W, 17V

Cascade Light Shield and Bracket Assembly

- 28.0503.01 Shield only
- 28.0503.02 Shield and nosepiece
- 28.1012.00 Shield bracket assy (2), Surf 4
- 28.1012.01 Shield bracket assy (2), White 2



- For A-dec 571 and Cascade lights

Pivot Stops

- 28.0679.01 (5)
for all lights after October 1987
- 28.0545.01 (5)
for all pre-Cascade unit- and post-mount lights
- 90.1044.00 Flexarm Stop Kit




28.0679.01



28.0545.01

Pre-Cascade, Cascade, and Performer On/Off Switch Kit


- 90.1039.00 Kit
- 28.0813.01 Long handle only, Surf 4
- 28.0813.02 Long handle only, Gray 3
- 28.1188.00 Short handle only



- Includes: toggle switch, (2) closed-end terminals and (2) On/Off switch handles

Pre-Cascade Dental Light Bulb


- 28.0682.00 Holder only
- 28.0704.00 Bulb and holder



- For pre-Cascade dental light
- Includes holder and halogen bulb 95W, 17V

Light Handle Kit


- 90.0367.01 Handle kit, Surf 4
- 90.0367.02 Handle kit, White 2



- Includes Loctite® and necessary hardware
- Includes both right and left handles

Light Switch Housing Replacement Kit


- 28.1464.01 Surf 4
- 28.1464.04 White 2



- Includes right housing, left housing, standoffs, screws, and pivot stop (28.0679.01)

Pre-Cascade and Cascade Lamp Socket Kit

- 90.0463.01



- To replace lamp socket for pre-Cascade and Cascade lights
- Includes lamp socket and wire assembly, and necessary hardware

A-dec 500 Light Post Trim Ring

- 77.0095.00 Surf 4
- 77.0095.01 White 2



Composants de l'éclairage 6300, transfo.

Cascade Light Cable Kit


90.1053.00 6300 light cable



- Replaces flexarm cable for all 6300 lights, prior to April 1, 2004
- Includes cable lube
- 137" long, cut desired length
- For Radius or Cascade lights

Circuit Breaker Kit Pre-Cascade and Cascade


28.0758.00



- Use to replace a 10A circuit breaker
- Includes circuit breaker, wire assemblies, and necessary hardware

Cascade and Pre-Cascade Intensity Switch

90.1043.00




- 6A, 24V
- For 6300 wall-, track-, Preference-, and ceiling-mount dental lights before serial number D415121
- For all Radius and Cascade lights

Pre-Cascade Dental Light Switch Service Kit

90.0372.00 Kit

28.0811.02 Cover only



- For all pre-Cascade lights except Performer
- Includes toggle switch, wire assembly, On/Off switch cover (Surf), and necessary hardware

Light Yoke Plug

28.1536.00 Surf 4


28.1536.01 White 2



- For A-dec 571 and 6300 lights

Transformer


90.1040.00



- For Cascade post, track, and ceiling lights before serial number D415121
- Pre-Cascade other than post-mount
- Replaces 28.0720.00
- 120V or 240V

Transformer


90.1041.00



- For Cascade Preference-mount lights and Cascade wall-mount lights before serial number D415121
- 120V or 240V
- Replaces 28.0724.00
- Includes hardware

Transformer

90.1042.00



- For Pre-Cascade post mount lights
- Includes transformer and hardware
- 120V or 240V

Câbles, transfo., interrupteurs et barres

Transformer

28.1588.00



- For all 6300 dental lights after serial number D415121 (2004)

Dental Light Circuit Board

28.1577.00 Kit



- For 6300 wall-, track-, Preference-, and ceiling-mount dental lights after serial number D415121 (2004)
- Includes bracket and standoffs

Intensity Switch

43.0054.00



- For 6300 wall-, track-, Preference-, and ceiling-mount dental lights after serial number D415121 (2004)
- A-dec 571 light

Lubriplate

048.107.00




- Use on light hub
- 3/8 ounce tube

Track Cap

028.008.00




CABLE ASSEMBLIES	PART NUMBER	PRICE
Track light track and trolley, after April 1, 2004	28.1604.00	
Track light flexarm, after April 1, 2004	28.1584.00	
Wall and Preference mount, after April 1, 2004	28.1585.00	
Post mount, after April 1, 2004	28.1586.00	
Ceiling mount, after April 1, 2004	28.1587.00	
Dental light 216" (5486 mm), A-dec 571	28.1838.00	

LIGHT POSTS		PART NUMBER WHITE	PART NUMBER SURF	PRICE
	for Performer: 24.5" (622 mm) H for installation on cuspidor or assistant's instrumentation: includes post, adapter, trim ring, and power cord	28.1262.01	28.1262.00	
	for Performer: 42.25" (1,073 mm) H x 17" (432 mm) W for installation on Performer chair; with chair only	28.1263.01	28.1263.00	
	for Pre-Cascade chairs: 2" (51 mm) OD x 40" (1,016 mm) H for use with chair-mount to mount A-dec post-mount lights to pre-Cascade chairs	28.1976.00	-	
		-	28.1143.00	
	for 1040 Chair: 2" (51 mm) OD x 46" (1,168 mm) H x 25.25" (641 mm) W only needed for installing A-dec post-mount light to Cascade 1040 chair	-	28.1105.00	
	for 1021/1221 Chair: 2" (51 mm) OD x 46" (1,168 mm) H x 25.25" (641 mm) W only needed for installing A-dec post-mount light to Cascade 1040 chair	-	28.0881.01	
	for Cascade post box (extension kit): 2" (51 mm) OD x 32" (813 mm) H for post-mount light to Cascade post box: includes power cord, post, and hardware	28.1030.01	28.1030.00	
	for A-dec 361 support center: 2" (51 mm) OD x 32" (813 mm) H for installation on A-dec 361 support center	99.1136.00	-	

Composants pour Performer

Toggle Switch Kit


90.1039.00



- 125V, SPST 20A
- Use small Surf toggle for Performer

Performer On/Off Switch Handle

28.1188.00



Washers for Performer

28.1175.01 Thrust 1.362 OD (5)

016.053.00 Thrust .250 ID



Thrust washer (5)
1.362 OD



Thrust .250 ID

Bearing, Compression Bolt, and Compression Spring


016.054.00 Bearing

28.1172.00 Bolt

013.100.00 Spring



Bearing, Thrust
.250 ID



Compression Bolt



Compression Spring
.480 OD x .500 LG

Performer Light Intensity Switch Kit

90.1045.00 Kit




- Includes switch and cable assembly

Performer Halogen Bulb, Socket, and Insulator Assembly

041.513.00 Bulb only

28.1289.00 Assembly



- Halogen bulb 55 W, 12V
- Assembly includes bulb, socket, and insulator
- Bulb available from your local automotive supplier, Type H3


Performer Bulb Cap

28.1213.00



Performer Light Shield

28.1166.00



Performer Light Cable Kit


90.1054.00 Cable kit



- Includes cable lube
- Replaces cable from light through flexarm to postbox

Protège-pieds

TOEBOARD COVERS BY MODEL	PART NUMBER	PRICE
A-dec 511		
511, seamless	62.0156.00	
511, sewn	62.0184.00	
A-dec 411		
411, seamless	62.0848.00	
411. sewn	62.0847.00	
A-dec 311		
311, seamless	62.0676.00	
311. sewn	62.0685.00	
Cascade		
1040, seamless	61.1679.00	
1040, sewn	61.2619.00	
Decade		
1011/1021/1221, seamless or sewn	61.1679.00	
1010 Square Toe	61.1421.00	
1010/1015 Thick	61.1593.00	
1020 Tubular	61.1632.00	
Performer		
I and II	61.1679.00	
III	61.2619.00	
Priority		
1005 Standard	61.0765.00	
EZ-Clean	61.0791.00	
Plus	61.0775.00	
EZ-Clean short	61.0797.00	
1005 Short	61.0765.02	




Selleries vis et molettes

DECADE AND CASCADE UPHOLSTERY SCREWS BY MODEL	PART NUMBER	PRICE
Decade		
Armrest (requires 3)	005.123.01	
Seat (requires 4)	61.2606.00	
Headrest	001.157.00	
Cascade		
Armrest	001.184.00	
Seat (requires 4)	61.2606.00	
Headrest	001.157.00	


Performer Upholstery Kits and Screws

005.155.00 Seat, each

61.2767.00 Backrest, each



- For mounting seat upholstery to frame
- Requires 4



- For mounting narrow or contoured backrest upholstery to frame
- Requires 4

A-dec 300 Double Articulating headrest knob Kit

90.1221.00



- Includes necessary hardware

A-dec 511 Glide Bar bracket kit


99.0719.00



- Includes necessary hardware

Sliding Wedge for Glide Bar Tension

61.1569.00




- Wear pad wedge for Cascade 1040
Decade 1221, 1021, 1020
Decade 1015/1010
Performer III

Headrest Locking Knob

61.1232.01 Surf 4

61.1232.02 Gray 3




- For Cascade, Decade, and Performer III double-articulating headrest (61.2265.00)
- For Decade single-articulating headrest (61.1227.00 and 61.1548.00)

Articulating Headrest Handle

61.1553.00 Handle

61.1551.00 Locking nut



- For double-articulating headrest
- Surf
- Replacing both handle and locking nut recommended

Remplacement des selleries

UN KIT SELLERIE COMPREND :

la tête, le dossier, l'assise, les accoudoirs et le protège pieds

UTILISEZ CE DOCUMENT COMME BON DE COMMANDE À RETOURNER À EUROTEC:

jean-philippe@eurotec-dental.fr

NOTE : Teintier page suivante

CHOIX DES COULEURS	PRODUCT DESCRIPTION	SEWN	SEAMLESS
Sellerie lisse			
☐ Amethyst	☐ Indigo	☐ Sky Blue	
☐ Black	☐ Lemongrass	☐ Tawny	
☐ Driftwood	☐ Paprika		
☐ Hazelnut	☐ Sapphire		
Sellerie avec coutures			
☐ Apricot	☐ Ebony	☐ Poppy	
☐ Arctic	☐ Fuchsia	☐ Riviera	
☐ Baltic	☐ Granite	☐ Sable	
☐ Campfire	☐ Lilac	☐ Sassy	
☐ Cave	☐ Lemon	☐ Sea Mist	
☐ Chamois	☐ Pacific	☐ Silver Pearl	
☐ Charcoal	☐ Papyrus	☐ Sky	
☐ Curry	☐ Parrot	☐ Sorbet	
☐ Cyan	☐ Pecan	☐ Timberwolf	
☐ Diplomat Blue	☐ Plum	☐ Vintage	
A-dec 511			
☐ White			
☐ Surf			
A-dec 311			
☐ Double Articulating Headrest			
A-dec 411			
☐ Knob-Style Headrest			
Performer III			
☐ Contoured Back			
Cascade/Decade/Decade Plus			
☐ Cascade 1040			
☐ Decade 1010			
☐ Decade 1011			

Remplacement des selleries

UPHOLSTERY COLOR CHOICES			PRODUCT DESCRIPTION		SEWN	SEAMLESS
Lisse ☐ Amethyst ☐ Indigo ☐ Sky Blue ☐ Black ☐ Lemongrass ☐ Tawny ☐ Driftwood ☐ Paprika ☐ Hazelnut ☐ Sapphire			Doctor's Stool ☐ 1601 Contoured			
			Assistant's Stool ☐ 1621 Contoured			
			☐ 1622 Contoured Round			
			☐ 1626 Contoured			
			PRODUCT DESCRIPTION		PART NUMBER	PRICE
Avec coutures ☐ Apricot ☐ Ebony ☐ Poppy ☐ Arctic ☐ Fuchsia ☐ Riviera ☐ Baltic ☐ Granite ☐ Sable ☐ Campfire ☐ Lilac ☐ Sassy ☐ Cave ☐ Lemon ☐ Sea Mist ☐ Chamois ☐ Pacific ☐ Silver Pearl ☐ Charcoal ☐ Papyrus ☐ Sky ☐ Curry ☐ Parrot ☐ Sorbet ☐ Cyan ☐ Pecan ☐ Timberwolf ☐ Diplomat Blue ☐ Plum ☐ Vintage			A-dec 311 Neck Support ☐ Cushion Only			
			Headrest Mechanisms ☐ Lever-Style (White)			
			☐ Lever-Style (Surf)			
			☐ Knob-Style			
			☐ Cascade/Decade/Performer/A-dec 200 (Surf)			
			☐ Cascade/Decade/Performer/A-dec 200 (White)			
			Performer III ☐ Slings*			
*Slings on Performer III chairs manufactured before April 2000 include integrated armrest caps. When ordering replacements for these slings, please order (2) P/N 61.2761.XX armrest caps.						

Remplacement des selleries

Teintes uniquement disponibles en option Luxe "Avec coutures"



cyan



sky




riviera




charcoal



lilac



granite



silver pearl




pecan



sea mist



chamois




curry



cave



lemon



arctic



fuchsia



sorbet




sassy



apricot



poppy



vintage


Teintes uniquement disponibles en sellerie standard "Lisse" et en option Luxe "Avec coutures"



baltic (sapphire)



diplomat blue (indigo)



plum (amethyst)



ebony (black)




sable (tawny)



papyrus (hazelnut)



campfire (paprika)



timberwolf (driftwood)



pacific (sky blue)



parrot (lemongrass)

Les couleurs lisses sont entre parenthèses.

Vérins

Lift Cylinders



Cascade, Decade and Performer III Chairs



Priority Chairs (S/N M2.0196 to end)



Priority Chairs (S/N M1.0001 to M2.0195)



A-dec 511 Chair

LIFT CYLINDERS MODEL	PART NUMBER	PRICE
Cascade, Decade, and Performer III Chairs	61.1287.00	
Priority Chairs (S/N M2.0196 to end)	61.0850.00	
A-dec 511 cylinder kit, Surf 4	90.1083.00	
A-dec 511 cylinder kit, White	90.1083.01	

Tilt Cylinders



Cascade Chairs



Performer III



Decade Chairs



A-dec 511 Chair


TILT CYLINDERS MODEL	PART NUMBER	PRICE
Cascade Chairs	61.2050.00	
Performer III	61.2050.01	
Decade Chairs	61.1267.00	
A-dec 511 Chair Kit, Surf	90.1085.00	
A-dec 511 Chair Kit, White	90.1085.01	

A-dec 300 Tilt Actuators




TILT ACTUATOR MODEL	PART NUMBER	PRICE
Actuator replacement kit	90.1228.00	

Performer I, II Lift, and Tilt Actuators



Lift Actuator



Tilt Actuator

LIFT ACTUATOR MODEL	PART NUMBER	PRICE
230V	61.2470.00	

TILT ACTUATOR MODEL	PART NUMBER	PRICE
230V	61.2410.00	

A-dec Manufactured Cylinders

À partir de la matière première, A-dec fabrique ses vérins. Une fois fabriqués, les vérins sont assemblés dans une pièce climatisée et hermétique afin d'éviter tout corps étranger. Testés à 1500 psi soit 100 bars, les vérins sont contrôlés afin d'assurer un fonctionnement sans faille.

A-dec garantie, pour votre tranquillité, ses vérins hydrauliques de dossier et de l'assise pendant 10 ans.

Électrovannes et huile

A-dec 300 and 400 Chair Solenoid Assembly

90.1225.00 Solenoid



- 21.6 VDC

A-dec 511 chair hydraulic Solenoid Kits

90.1071.00 Red wires, 220-240V



Priority Chair Solenoids

60.0601.00 White, 24VDC



- 24 VDC
- O-ring packages:
For two-piece block (035.041.02)
For one-piece block (030.015.02)

Cascade, Decade, Performer III, and Positioner Chair Solenoids


61.1337.00 Red wires, 240V



- O-ring packages:
For two-piece block (035.041.02)
For one-piece block (030.015.02)

Hydraulic Oil

61.0197.00




- Hydraulic oil to add to all A-dec hydraulic chairs
- Includes noise eliminating additive
- 1 pint (16 ounces)

Flow Adjustment Screw

61.0460.00 Screw with O-ring

030.004.02 O-ring (10)



- For hydraulic manifold on all chairs
- Flow adjustment screw with O-ring

Truss Head Slot Screw Button Head Socket Screw

001.002.01 Truss (5)

002.118.02 Button (5)



Truss



Button

- Use on all two piece hydraulic manifolds to retain flow adjustment screws
- #2-56 x 1/8
- Use on all one piece hydraulic manifolds to retain flow adjustment screws
- 4-40 x 3/16

Condensateurs, alimentations et platines

Printed Circuit Board Kit	PART NUMBER	PRICE
A-dec 511 Chair	90.1072.00	

Circuit Board	PART NUMBER	PRICE
Cascade, Decade, and Performer II, III Replaces 90.1029.00, 90.1029.01, 61.1214.01, 61.1373.01 and 61.2510.00, 61.1217.01 and 61.2512.00	90.1029.02	


Chair Power Supplies	PART NUMBER	PRICE
A-dec 511 Chair		
300W, 220-240VAC	90.1075.00	

Remote Master Toggle to Power Supply Kit
90.1129.00



- Includes 25' (7.62 m) of tubing, QD and Tee connectors
- For connecting master toggle to the chair power supply

Capacitors for Cascade, Decade, Positioner, Priority, Performer I, II, and III Chairs




DESCRIPTION	PART NUMBER	PRICE
2.0 Dia Boot	041.529.00	
15 UF Capacitor	041.504.00	
5 - 7 1/2 UF capacitor with boot for 250v chairs	90.1036.00	

Replacement Capacitors	PART NUMBER	PRICE
A-dec 311 and 411 Capacitor Kits		
Replacement Kit, 220-240v	90.1200.00	
A-dec 511 Capacitor		
60 UF Capacitor, 200v	041.642.00	
50 UF Capacitor, 220v	041.643.00	
15 UF Capacitor, 420v	041.644.00	

Capteurs de position

Limit Switch


044.184.01



- Base up function for:
Cascade 1040
Decade 1020/1021
Decade 1010/1011
Positioner 1030
Priority 1005
Performer III
- Back up function for:
Positioner 1030
Priority 1005
- Chair stop function for:
Decade 1020/1021
Positioner 1030
Priority 1005

Limit Switch


61.2065.00



- Back up function for:
Cascade 1040

Limit Switch


044.183.00



- Chair stop function for:
Cascade 1040,
Decade 1021 after June 1998,
and Performer I and II

Limit Switch


044.184.00



- Base up and base down
function for: Performer I and II
- Back up function for:
Decade 1020/1021
Decade 1010/1011
Performer III

Limit Switch

61.2415.00



- Back up function for:
Performer I only

A-dec 500 Positioning Sensor Kit

90.1069.00 Position sensor kit




- For base or back positioning
A-dec 500 chairs

Stop Plate Switch Replacement Kit

90.1095.00 Surf 4

90.1095.01 White 2



- For use with A-dec 511


A-dec 300/400 Positioning Sensor Kit

90.1183.00 Position sensor kit

Potentiomètre et butée mécanique

Potentiometer

041.372.00




- For Cascade, Decade, and Performer II and III Chairs

Base Down Stop for Chairs*

61.2576.00 Radius light

61.3053.00 Radius monitor mount




- For Radius light 1³/₈" is anodized aluminum
- For Radius monitor mount 1¹/₂" is black anodized aluminum

NOTE: *You must install the appropriate base down stop when you install a Radius dental light or monitor mount on a chair.
Omitting the base down stop causes damage to the motor-pump cover when the chair is moved to full base down.
Reprogram the chair entry/exit function after installing the stop.

Freins de rotation, câbles électriques et cordons d'alimentation


Cascade Chair Swivel Brake Kit

90.1016.00 Swivel brake kit



Decade Chair Swivel Brake Kit


61.1537.01 Pads (2)



Decade Swivel Brake Kit

61.3064.00 Assembly, Surf

61.1537.01 Pads (2)



Priority Chair Brake Cable

Brake cable

PART NUMBER

61.0025.00

PRICE

Brake cable mount

61.0081.00

Cascade/Decade Electrical Wiring/Cables

PART NUMBER

PRICE

Back up, base and safety plate limit switch for 1015/1020/1021/1040

61.2099.00

Base electrical wiring cable for potentiometer

61.1502.00

Back electrical wiring cable for potentiometer

61.1503.00

Priority Electrical Wiring/Cables

PART NUMBER

PRICE

Main cable, M1-0001 through D3-0134

61.0612.00

Main cable, D3-0135 through H3-0067

61.0599.00

Main cable, H3-0038 through present

61.0668.00

Stop plate cable, S/N M1-0001 through present

61.0614.00

Accessory Power Cord

43.0036.00 Power cord without plug 511, 240v

28.1276.00 Power cord without plug for Decade, Cascade, and Performer, 230v




Pédales électriques

8-Button Footswitch

62.0163.01 Programmable Footswitch, Gray 3

For use on all current A-dec chairs




8-Function Footswitch*

Service Parts Only

8 Function Footswitch (Part No. 61.2027.00) replaced by 8-button Footswitch in 1997

*Base and Cover Assy Not Available For Sale



Barres de montage et accoudoirs

2" Post Adapter to 511 Chair

77.0921.00



- Mounts an A-dec 300 2" post-mount delivery and support center to an A-dec 511 chair
- Includes necessary hardware

2" Post Adapter to 411 Chair

77.1301.00



- Mounts an A-dec 300 2" post-mount delivery and support center to an A-dec 411 chair
- Includes necessary hardware

2" Post Adapter to Performer Chair

77.0922.00



- Mounts an A-dec 300 2" post-mount delivery and support center to a Performer III 8000 chair
- Includes necessary hardware

2" Post Adapter to Cascade 1040 Chair

35.1742.00



- Mounts an A-dec 300 2" post-mount delivery and support center to a Cascade 1040
- Includes necessary hardware

2" Post Adapter to Priority 1005 Chair

35.1743.00



- Mounts an A-dec 300 2" post-mount delivery and support center to a Priority 1005 chair
- Includes necessary hardware

A-dec 500 Quick disconnect Kit

90.1138.00 Chair mount QD kit



- For A-dec 511 chair-mounting or A-dec 561 support center-mounting

2" Post Adapter to Decade 1021/1221 Chair

35.1741.00



- Mounts an A-dec 300 2" post-mount delivery and support center to a Decade 1011/1021/1221 chair
- Includes necessary hardware

A-dec 311 1/2" Hole Plug

018.002.06 White 2



Performer Pivot Armrest

61.3125.00




Barres de montage et accoudoirs

Cascade or Radius Delivery System to Cascade 1040 Chair Mount

35.1372.00 Cascade/Cascade

35.1572.00 Radius/Cascade



- Includes necessary hardware

Decade Chair Mount

35.1150.00 1010/1015/1020

35.1150.01 1011/1021/1221



- Mounts pre-Cascade 2" post to Decade chairs
- Includes necessary hardware

Cascade 1040 Chair Adapter


35.1465.00



- Adapts Decade 2" post to 1040 chair

Radius Delivery to Decade 1021 Chair Mount


35.1667.00 Radius/Decade



- Includes necessary hardware
- For external umbilical

Cascade Delivery System to Decade 1021 Chair Mount


35.1401.00



- Includes necessary hardware

Priority Chair Mount

35.0824.02



- Mounts pre-Cascade 2" post to Priority 1005 chair
- Includes necessary hardware

Cascade Adapter Insert

35.1405.00




- For use with chair mount to adapt Cascade unit to non A-dec chair

Molettes et roulettes avant 2015

Cart Height Adjustment Knob

027.007.02 Dark Surf



- Knob to adjust height on "H-frame" adjustable height carts

Stool Back Height Adjustment Knob


027.063.00 1622 stools, Surf 4

027.063.01 1622 stools, White 2

027.065.00 1601 stools, Surf 4

For stools manufactured before March 24, 2003

027.024.01 1600 stools, Gray



- Knob to adjust height on the back of the stool or the torso bar

Stool Back Height Adjustment Knob

027.074.00 Surf 4


027.074.01 White 2



- Use for 1601 stools manufactured after March 24, 2003
- With molded-in bolt

Stool Height Adjustment Knob

027.035.01 Gray



- Knob to adjust height on the assistant's stool with rotating torso


One-Piece Casters for Carts and Stools

65.1625.03 Tile (5), Surf

65.1625.04 Tile (5), Gray 3

65.0286.02 Carpet (5), Surf

65.0286.05 Carpet (5), Gray 3



- Tile caster has plastic wheel with rubber surface for use on tile floors
- Carpet caster has all plastic wheel surface for use on carpeted floors


Caster for Oregon Cart

042.615.02 Tile, Gray 3



Caster for U-frame Carts


16.0081.00 each



- 2" ball caster for A-dec U-frame carts
- Includes screw and washer

Caster for Carts (beginning June 2014)

30.0861.00 Tile, Surf



Vérins et dossiers avant 2015

Cylinders for Doctor's Stools



Cascade Doctor's Stools, Short
Height range: 15.25" to 18.5"



Cascade Doctor's Stools, Standard
Height range: 18.5" to 25.0"

CYLINDER MODEL	PART NUMBER	PRICE
Cascade Doctor's Stools, Short Height range: 15.25" to 18.5"	65.0227.00	
Cascade Doctor's Stools, Standard Height range: 18.5" to 25.0"	65.0208.00	

Cylinders for Assistant's Stools




Cascade Assistant's Stools Pre-Cascade Assistant's Stools
(Series 1620 from 06/93 through present)

CYLINDER MODEL	PART NUMBER	PRICE
Cascade Assistant's Stools Pre-Cascade Assistant's Stools (Series 1620 from 06/93 through present)	65.0209.00	

Contour Stool Back Cover Kit

65.0289.00 Surf 4

65.0289.01 White 2



- Includes back cover, trim cover, and screws

Supports plateaux

A-dec 532/542 Traditional Tray Holders

- Arm-mounted rotating tray holder

PRODUCT DESCRIPTION	PART NUMBER SURF	PART NUMBER WHITE	PRICE
Standard tray holder 10.3" x 14.3" (262 mm x 363 mm)	-	77.1024.00	
Large tray holder 12.5" x 17.3" (318 mm x 439 mm)	77.0295.00	77.1023.00	



Traditional Standard
Tray Holder (for A-dec 532/542)



Traditional Large
Tray Holder (for A-dec 532/542)

A-dec 541/545 Tray Holders

- Arm-mounted rotating tray holder

PRODUCT DESCRIPTION	PART NUMBER SURF	PART NUMBER WHITE	PRICE
Standard tray holder 9.35" x 12.66" (278 mm x 322 mm)	77.0462.00	77.0462.01	
Large tray holder 12.5" x 17.3" (318 mm x 439 mm)	77.0543.00	77.0543.01	



Standard Tray Holder
(for A-dec 541/545)



Large Tray Holder
(for A-dec 541/545)



Standard Tray Holder
(for A-dec 581)



Large Tray Holder
(for A-dec 581)

Supports plateaux

A-dec 533 Continental Tray Holders

- Arm-mounted rotating tray holder

PRODUCT DESCRIPTION	PART NUMBER SURF	PART NUMBER WHITE	PRICE
Standard tray holder 10.3" x 14.3" (262 mm x 363 mm)	-	77.1024.01	
Large tray holder 12.5" x 17.3" (318 mm x 439 mm)	77.0295.01	77.1023.01	



Continental Standard
Tray Holder (for A-dec 533)



Continental Large
Tray Holder (for A-dec 533)

A-dec 332/334 Traditional Tray Holders

- Arm-mounted rotating tray holder

PRODUCT DESCRIPTION	PART NUMBER	PRICE
Standard tray holder 10.3" x 14.3" (262 mm x 363 mm)	77.0933.00	
Large tray holder 12.5" x 17.3" (318 mm x 439 mm)	77.0932.00	



Standard tray holder
(for A-dec 334)



Large tray holder
(for A-dec 334)

A-dec 333/335 Continental tray holders

- Arm-mounted rotating tray holder

PRODUCT DESCRIPTION	PART NUMBER	PRICE
Standard tray holder 10.3" x 14.3" (262 mm x 363 mm)	77.0942.00	
Large tray holder 12.5" x 17.3" (318 mm x 439 mm)	77.0943.00	



Standard tray holder
(for A-dec 335)



Large tray holder
(for A-dec 335)

Plateaux et tapis


Stainless Steel Trays	PART NUMBER	PRICE
Stainless steel tray, standard 9.5" x 13.5" (241 mm x 343 mm)	043.001.00	
Stainless steel tray, large 11.5" x 17" (292 mm x 432 mm)	043.006.00	



Washable Tray Pad

- Firmly holds stainless instrument tray to flat surface
- Large size: 6" x 11.5" (152 mm x 292 mm)
- Small size: 6" x 9" (152 mm x 229 mm)

PRODUCT DESCRIPTION	PART NUMBER	PRICE
Washable tray pad, large (5)	39.0790.02	
Washable tray pad, small (5)	39.0790.03	



Performer Tray Holder

- Includes tray holder, mounting bracket, and mounting hardware

PRODUCT DESCRIPTION	PART NUMBER	PRICE
Standard tray holder 9.5" x 13.5" (241 mm x 343 mm)	39.1380.00	
Large tray holder	39.1390.00	
Thumb screw, white	027.070.00	




- Shown on control head with optional tray


Tray Holders

- Arm-mounted rotating tray holder


PRODUCT DESCRIPTION	PART NUMBER	PRICE
For Radius 2132, Cascade 2150, 2151, and 2152 10" x 14" (254 mm x 356 mm) (includes tray and pad) Tray Color: White	39.1115.00	
For Radius 2122, Cascade 2140, 2141, and 2142 (serial numbers E443922 and subsequent) 10" x 14" (254 mm x 356 mm) Tray Color: White	17.1051.00	
For Cascade 3171, 3175, 4631, 4635, 2615, 2671, 318 19.5" x 12.5" (241 mm x 318 mm)	17.1047.00	



39.1115.00



17.1047.00




17.1051.00

12 O'Clock Worksurface Bumper and Plug

77.0541.01 Bumper, white


77.0542.01 Plug, white



- For worksurface when no tray is present

Radius Flat-Panel Mount

35.1723.00



- Requires Radius monitor post 35.1722.00 or Performer monitor post 83.0486.00

Performer Monitor Post Assembly

83.0486.00




- For use with Radius Flat-Panel Mount, 35.1723.00

Protections units

Instrument Pads

75.0085.01 Cascade 2132 (100)



- Continental handpiece holder
- Package of 100 plastic disposable covers

77.0070.00 A-dec 533, Surf 4

77.0070.01 A-dec 533, White 2



- Easily removable
- Dishwasher safe


77.1225.00 A-dec 333/335, White 2



- Easily removable
- Dishwasher safe

Cascade & Radius Brake Handle Kit


90.1013.00 Brake handle



- Includes brake handle and necessary hardware
- Monts for right- or left-hand use

Control Top Assembly

39.1024.00



Index des références

Index	Page(s)	11.1025.02	33	12.0910.06	33	12.1258.00	32, 33, 68	17.0594.02	72
0.01862301	1	11.1025.12	33	12.0910.12	33	013.003.01	44, 45	17.0600.01	72
001.002.01	44, 91	11.1027.00	31, 33, 68	12.0913.00	39, 74, 75	013.004.00	76	17.1047.00	103
001.157.00	86	11.1027.01	31, 68	12.0915.00	74, 75	013.015.00	55	17.1051.00	103
001.184.00	86	11.1074.00	31, 33, 68	12.0920.00	38	013.025.00	56	018.002.02	20
001.191.00	60	11.1127.01	33	12.0931.00	78	013.027.00	55	018.002.04	12, 20
001.253.00	5	11.1131.00	33	12.0934.00	39, 74, 75	013.064.01	44, 45, 46	018.002.05	12, 20
002.032.00	96	11.1178.00	33	12.0953.00	74, 75	013.072.00	54	018.002.06	12, 20, 97
002.118.02	91	11.1191.00	72	12.0954.01	38	013.100.00	84	018.002.08	12, 20
002.120.00	95	11.1192.00	72	12.0977.01	76	013.123.00	59	018.010.02	20
002.135.00	5	11.1194.00	72	12.0982.01	76	13.0350.00	38	018.010.03	20
003.011.01	20	11.1235.01	32, 33, 68	12.0983.00	76	13.0361.00	38	018.026.02	20
003.071.00	5	11.1272.00	31	12.0985.00	76	13.0375.00	17	018.026.03	20
003.078.00	52	11.1280.01	70, 71	12.0986.00	76	13.0381.00	17	018.026.04	20
004.005.01	16, 57	11.1280.02	70, 71	12.0987.00	76	13.0384.00	17	018.026.05	20
004.005.02	7, 15, 16, 57	11.1283.01	32	12.0988.00	76	014.006.00	20	018.086.00	76
004.008.01	15	11.1283.03	32	12.0989.00	76	014.007.00	20	018.094.00	6
004.073.01	96	11.1284.01	70, 71	12.0991.00	76	14.0308.02	7	018.119.00	12
004.137.00	7	11.1286.00	31	12.1013.00	76	14.0326.00	57	19.0001.03	57
005.012.03	96	11.1286.01	31	12.1024.00	74, 75	14.0327.01	7	19.0100.00	21
005.110.00	4	11.1287.00	34	12.1028.00	74	14.0332.01	7	19.0220.01	21
005.123.01	86	11.1287.01	34	12.1031.00	74	14.0408.00	7	021.016.04	16
005.155.00	86	11.1288.00	31	12.1034.00	74	14.0416.00	7, 76	021.027.01	15
006.009.00	23, 36, 37, 38, 40, 42, 78	11.1288.01	31	12.1054.02	74	14.0424.00	7	021.034.01	16
006.069.00	36, 37, 38, 40, 42, 78	11.1292.00	31	12.1070.01	37	14.0453.00	4	022.009.01	16
006.133.01	74	11.1292.01	31	12.1071.02	37	14.0455.02	5	22.0060.00	50
007.001.00	42	11.1293.00	34	12.1087.00	33	14.0455.04	5	22.0110.00	52
008.023.00	24	11.1293.01	34	12.1088.00	33, 68	14.0457.00	4	22.0120.00	52
009.005.00	24	11.1294.00	31	12.1093.00	32, 33, 68	14.0458.00	4	22.0440.01	8
009.018.00	24	11.1294.01	31	12.1109.01	34, 68	14.0458.01	4	22.0440.02	8
009.037.00	24	11.1296.00	31	12.1116.00	34, 68	14.0464.00	5	22.0462.02	42
010.002.00	4	11.1296.01	31	12.1121.00	31, 34	14.0464.02	5	22.0462.03	42
10.0014.01	69	11.1297.00	34	12.1123.00	72	14.0466.00	4	023.001.02	15
10.0015.00	69	11.1300.01	31	12.1132.00	34	14.0466.01	4	023.001.03	7, 15
010.045.02	76	11.1300.03	31	12.1228.00	32	14.0467.00	4	023.004.03	15
010.073.01	48	11.1301.03	31, 34, 68	12.1228.01	32	14.0468.00	4, 5	023.004.04	15
011.038.01	42, 44, 45	11.1302.01	32	12.1232.00	32	14.0468.01	4, 5	023.014.02	15
11.0070.00	69	11.1302.03	32	12.1232.01	32	14.0471.00	5	23.0018.00	47
11.0090.00	69	11.1305.00	31	12.1236.00	32	14.0484.00	5	23.0019.00	47
11.0170.00	69	11.1305.01	31	12.1236.01	32	14.0491.01	6	023.019.02	14
11.0440.00	69	11.1307.00	32	12.1237.00	31	016.044.00	18	023.019.03	14
11.0490.00	69	11.1307.01	32	12.1237.01	31	016.053.00	84	23.0020.02	47
11.0998.01	31, 33, 68	11.1315.00	34	12.1238.00	34	016.054.00	84	023.025.00	14
11.1007.00	72	11.1335.00	31, 33, 68	12.1238.01	34	16.0081.00	99	23.0027.00	47
11.1014.12	33	12.0677.01	77	12.1239.00	31	016.108.00	18	023.028.02	15
11.1016.00	72	12.0709.00	78	12.1239.01	31	017.007.01	52, 96	023.029.00	14
		12.0783.00	77	12.1246.00	33	17.0593.01	72	023.030.02	16
		12.0784.00	77	12.1257.00	34, 68	17.0594.01	72	023.031.01	16

Index des références

023.036.00	15, 57	24.0118.00	23	024.202.01	35, 52	025.109.00	16	28.0881.01	83
023.045.04	14	24.0128.01	22	024.202.02	35, 52	025.109.01	16	28.1004.00	81
023.057.01	16	24.0129.01	22	024.203.01	18	026.021.01	21	28.1004.01	81
023.065.01	15	24.0136.01	8	024.203.02	35	026.036.00	42	28.1012.00	81
023.070.00	15, 57	24.0137.01	8, 39	024.209.01	35	026.062.00	23, 40	28.1012.01	81
023.079.00	14	24.0138.01	22	024.209.02	35	026.066.00	21	28.1030.00	83
023.088.00	29	024.147.00	30	024.210.01	35	026.073.02	41	28.1030.01	83
023.094.01	16	024.147.01	30	024.210.02	35	026.074.01	41	28.1105.00	83
023.095.01	16	024.148.01	30	024.214.01	35, 52	026.118.00	23, 24	28.1143.00	83
23.0173.00	42	024.148.02	30	024.214.02	35, 52	026.121.00	17	28.1166.00	84
23.0220.02	47	024.148.04	30	024.219.01	35	026.123.00	37	28.1172.00	84
023.804.00	15	024.151.01	35, 52	024.219.02	35	026.125.00	17	28.1175.01	84
023.805.01	15	024.151.03	35, 52	024.220.01	35	026.130.00	17	28.1188.00	81, 84
023.814.01	14	024.152.01	35	024.222.01	35	026.132.00	17	28.1213.00	84
23.0872.01	43	024.152.03	35	024.222.02	35	026.134.00	17	28.1262.00	83
23.0872.02	43	024.153.18	35	024.223.02	48	026.157.00	17	28.1262.01	83
23.1011.00	44	024.153.20	35	024.223.12	48	026.158.00	23, 40	28.1263.00	83
23.1015.00	48	024.155.02	48	024.228.01	18	026.160.01	53	28.1263.01	83
23.1021.01	44	024.155.12	48	024.228.02	35	026.161.00	29	28.1276.00	95
23.1064.01	44, 45	024.156.09	35	024.229.01	35, 52	026.174.01	28	28.1289.00	84
23.1067.00	44, 45	024.156.11	35	024.229.02	35, 52	026.175.01	28	28.1436.00	64
23.1068.01	44, 45	024.162.01	35	024.232.01	35	026.186.00	21	28.1464.01	81
23.1090.00	43, 44, 45	024.162.03	35	24.0234.01	23	27.0005.03	41	28.1464.04	81
23.1096.00	45	024.166.01	35	024.244.02	35	027.007.02	99	28.1481.00	64
23.1101.01	43	024.166.03	35	024.244.03	35	027.008.00	37, 42	28.1501.00	11
23.1110.00	44, 45	024.176.06	35	24.0251.01	23	027.017.02	56	28.1501.01	11
23.1112.00	43, 44, 45	024.176.08	35	24.0363.00	22	027.018.02	56	28.1536.00	82
23.1150.00	44, 45	024.177.01	35	24.0363.04	22	027.024.01	99	28.1536.01	82
23.1193.01	44, 45	024.177.03	35	24.0382.00	23	027.035.01	99	28.1577.00	63, 83
23.1194.00	44, 46	024.179.01	35, 52	24.0388.02	7	027.057.00	58	28.1584.00	83
23.1195.00	44	024.179.03	35, 52	24.0410.00	60	027.063.00	99	28.1585.00	83
23.1208.00	48	24.0182.02	23	24.0469.00	22	027.063.01	99	28.1586.00	83
23.1208.01	48	024.184.01	35	24.0470.00	22	027.065.00	99	28.1587.00	83
23.1232.01	44, 45, 46	024.184.03	35	24.0471.00	22	027.070.00	103	28.1588.00	83
23.1234.00	46	024.190.01	35	24.0474.00	18, 19	027.074.00	99	28.1604.00	83
23.1235.01	46	024.190.02	35	24.0475.00	18, 19	028.008.00	72, 83	28.1624.00	80
23.1247.00	45	024.191.05	35	025.002.01	24	28.0503.01	81	28.1634.00	80
23.1302.00	44	024.191.07	35	025.003.01	24	28.0503.02	81	28.1635.00	80
23.1304.00	45	024.192.01	35	025.007.01	16	28.0545.01	81	28.1832.00	24
23.1305.00	45	024.192.03	35	025.007.02	16	28.0679.01	81	28.1862.00	79
23.1309.00	44	024.193.01	35	025.015.00	16	28.0682.00	81	28.1924.00	79
024.015.04	26	024.193.03	35	025.015.01	16	28.0704.00	81	28.1925.00	79
024.015.05	26	024.197.01	35	025.015.02	16	28.0720.00	82	28.1976.00	83
024.090.04	27	024.197.02	35	025.017.01	24	28.0758.00	82	29.0050.01	8
024.090.05	27	024.198.01	35	025.019.01	24	28.0811.02	82	29.0100.00	57
24.0100.03	42	024.198.02	35	025.036.01	24	28.0813.01	81	29.0146.00	41
24.0105.01	57	024.200.00	5	025.052.00	42	28.0813.02	81	29.0185.00	41

Index des références

29.0830.00.....	42	33.0101.02.....	57	35.1150.01.....	98	036.052.04.....	27	38.0540.00.....	54
29.0831.00.....	56	33.0103.00.....	56	35.1372.00.....	98	036.053.03.....	27	38.0555.00.....	54
030.001.02.....	10	33.0104.05.....	40	35.1401.00.....	98	036.053.04.....	27	38.0604.00.....	51
030.002.02.....	10, 44, 45, 78	33.0108.02.....	36, 78	35.1405.00.....	98	036.054.03.....	27	38.0607.01.....	50
030.003.02.....	10, 54	33.0131.02.....	37	35.1465.00.....	98	036.054.04.....	27	38.0610.00.....	51
030.004.02.....	10, 56, 57, 58, 91	33.0132.00.....	41, 61	35.1572.00.....	98	036.066.03.....	27	38.0612.00.....	51
030.005.02.....	10	33.0134.00.....	41, 51, 61, 74, 75, 78	35.1667.00.....	98	036.066.04.....	27	38.0612.01.....	51
030.006.02.....	4	33.0138.00.....	41, 51, 61, 74, 75, 78	35.1722.00.....	103	036.103.03.....	27	38.0631.00.....	55
030.008.02.....	10	33.0160.00.....	74, 75	35.1723.00.....	103	036.103.04.....	27	38.0638.00.....	55
030.009.02.....	10	33.0168.00.....	36	35.1741.00.....	97	036.114.01.....	27	38.0648.00.....	55
030.011.02.....	10, 74, 75, 77	33.0168.01.....	36	35.1742.00.....	97	036.114.03.....	27	38.0659.00.....	50
030.012.02.....	10	33.0173.00.....	37	35.1743.00.....	97	036.116.01.....	27	38.0711.01.....	8, 58
030.013.02.....	10	034.001.01.....	9, 54	036.003.03.....	26	036.116.03.....	27	38.0712.00.....	58
030.014.02.....	10, 72	034.002.01.....	9	036.003.04.....	26	38.0054.01.....	8	38.0713.00.....	58
030.015.02.....	10, 91	034.003.01.....	9, 43	036.004.03.....	26	38.0054.02.....	2, 8, 56	38.0717.00.....	58
030.016.02.....	10, 73	034.004.01.....	10	036.004.04.....	26	38.0075.00.....	52	38.0718.00.....	58
030.019.03.....	10	034.012.01.....	9, 32, 33, 68	036.005.03.....	26	38.0075.03.....	52	38.0729.00.....	58
030.021.02.....	10	034.013.01.....	9, 31, 33, 68	036.005.04.....	26	38.0086.00.....	41, 78	38.0761.00.....	50
030.023.02.....	10	034.014.01.....	9, 31, 33, 68	036.006.03.....	26	38.0098.00.....	37	38.0763.00.....	51
030.026.01.....	10, 72	034.018.02.....	9, 34, 68	036.006.04.....	26	38.0178.00.....	56	38.0766.02.....	58
030.027.01.....	10, 72	034.019.01.....	9, 34, 68	036.009.04.....	26	38.0179.02.....	8, 56	38.1601.00.....	49
030.034.02.....	10, 72	34.0031.00.....	40	036.009.05.....	26	38.0181.00.....	56	38.1601.01.....	49
030.035.02.....	10, 70, 71	34.0032.00.....	40	036.010.03.....	26	38.0186.01.....	8	38.1603.00.....	49
030.224.01.....	10, 73	34.0033.00.....	40	036.010.04.....	26	38.0186.02.....	8	38.1700.00.....	49
030.228.02.....	10, 77	34.0036.00.....	40	036.011.03.....	26	38.0193.01.....	56	38.1764.00.....	49
30.0861.00.....	99	34.0037.00.....	40	036.013.03.....	26	38.0194.01.....	56	38.1775.00.....	58
031.130.01.....	4, 10	34.0037.01.....	40	036.013.04.....	26	38.0201.00.....	50	38.1780.00.....	53
033.003.01.....	9	34.0043.00.....	39	036.014.02.....	26	38.0223.02.....	56	38.1783.00.....	53
33.0008.02.....	36	034.107.01.....	9, 32, 33, 68	036.018.03.....	26	38.0234.02.....	56	38.1785.01.....	53
33.0009.04.....	36	034.213.00.....	5	036.018.04.....	26	38.0237.00.....	52	38.1787.01.....	8, 53
33.0025.00.....	41, 61	035.003.01.....	10	036.019.03.....	26	38.0316.00.....	57	38.1802.00.....	51
33.0025.01.....	41, 61, 78	035.004.01.....	10	036.019.04.....	26	38.0320.02.....	51	38.1802.01.....	51
33.0031.02.....	42	035.024.01.....	10	036.020.03.....	26	38.0321.02.....	51	38.1805.00.....	49
33.0034.02.....	36	035.025.01.....	10	036.031.02.....	27	38.0375.00.....	55	38.1805.01.....	49
33.0036.00.....	42	035.030.01.....	9, 43	036.031.04.....	27	38.0467.00.....	56	38.1808.01.....	53
33.0036.01.....	42	035.032.01.....	10	036.032.02.....	27	38.0504.07.....	55	38.1823.00.....	53
33.0036.02.....	42	035.033.01.....	9, 54	036.032.04.....	27	38.0504.08.....	55	38.1824.00.....	53
33.0036.03.....	42	035.034.01.....	2, 10, 54	036.033.02.....	27	38.0504.09.....	55	38.1825.00.....	53
33.0037.02.....	42	035.041.02.....	10, 91	036.033.04.....	27	38.0507.01.....	8	38.1829.01.....	49
33.0048.05.....	37	035.049.01.....	9	036.034.02.....	27	38.0510.00.....	2, 54	38.1832.00.....	51
33.0076.01.....	36	035.050.01.....	9	036.034.04.....	27	38.0516.00.....	53, 54	38.1834.00.....	49
33.0077.01.....	36	035.053.01.....	10, 70	036.043.01.....	26	38.0517.00.....	2, 54	38.1838.00.....	49
33.0080.03.....	37	035.062.01.....	43	036.043.02.....	26	38.0518.00.....	54	38.1862.00.....	49
33.0084.01.....	38	35.0574.01.....	8	036.044.03.....	26	38.0519.01.....	8	39.0790.02.....	103
33.0085.02.....	38	35.0824.02.....	98	036.044.04.....	26	38.0520.00.....	54	39.0790.03.....	103
33.0096.02.....	37	35.1150.00.....	98	036.051.03.....	27	38.0525.00.....	54	39.0820.00.....	57
33.0100.02.....	57			036.051.04.....	27	38.0526.00.....	53, 54	39.1050.00.....	55
				036.052.03.....	27	38.0537.01.....	55	39.1051.00.....	55

Index des références

39.1054.00	55	41.1454.01	18	052.015.00	7	61.2050.00	90	75.0035.03	73
39.1115.00	103	41.1455.00	19	052.023.00	60	61.2050.01	90	75.0039.00	74, 75
39.1380.00	103	41.1456.00	19	54.0162.03	1	61.2056.00	95	75.0078.00	72
39.1390.00	103	41.1457.00	18	54.0177.00	2	61.2065.00	93	75.0083.00	60
40.0145.00	68	41.1464.01	18	54.0451.01	1	61.2099.00	95	75.0084.00	81
40.0250.00	67	41.1467.00	18	54.0478.00	1	61.2189.00	96	75.0084.01	81
40.0849.00	77	41.1477.00	21	54.0480.01	1	61.2227.00	95	75.0085.01	104
41.0062.01	20	41.1490.00	19	54.0481.01	1	61.2228.00	95	75.0085.03	77
41.0074.00	20	41.1491.00	19	56.0339.00	42	61.2265.00	86	75.0091.00	74
41.0192.02	20	41.1506.00	18	060.001.00	25	61.2410.00	90	75.0145.01	74
41.0192.03	20	42.0008.00	77	060.002.00	25	61.2415.00	93	75.0150.00	12, 80
41.0217.00	21	42.0085.00	77	60.0601.00	91	61.2470.00	90	75.0911.01	8, 55
041.372.00	94	42.0090.00	77	61.0025.00	95	61.2510.00	92	75.0948.00	20
41.0407.00	19	42.0534.00	39	61.0081.00	95	61.2512.00	92	75.2018.00	12, 80
41.0408.00	19	42.0537.04	78	61.0197.00	91	61.2576.00	94	75.2251.00	80
41.0411.00	20	042.615.02	99	61.0460.00	91	61.2606.00	86	76.2011.00	76
41.0413.00	19	42.0683.04	36	61.0599.00	95	61.2619.00	85	76.2020.02	76
41.0413.01	19	42.0701.00	42	61.0612.00	95	61.2767.00	86	76.2025.00	17
41.0414.00	19	043.001.00	103	61.0614.00	95	61.3041.00	96	76.2030.00	17
41.0416.00	18	043.006.00	103	61.0668.00	95	61.3048.00	96	77.0021.01	11
41.0417.00	18	43.0010.00	73	61.0765.00	85	61.3049.00	96	77.0021.03	11
041.427.01	62	43.0012.00	83	61.0765.02	85	61.3053.00	94	77.0042.01	73
041.504.00	92	43.0036.00	95	61.0775.00	85	61.3064.00	95	77.0042.04	73
041.513.00	84	43.0054.00	62, 83	61.0791.00	85	61.3066.00	95	77.0043.01	73
041.529.00	92	43.0057.00	62	61.0797.00	85	61.3125.00	97	77.0043.04	73
041.554.00	63	43.0092.01	65	61.0850.00	90	62.0156.00	85	77.0044.00	73
041.555.00	63	43.0093.00	65	61.1214.01	92	62.0163.01	96	77.0044.01	73
041.582.00	19	43.0095.00	65	61.1217.01	92	62.0184.00	85	77.0063.00	59
041.637.00	63	43.0328.01	65	61.1227.00	86	62.0186.00	96	77.0063.01	59
041.637.02	63	43.0329.00	65	61.1232.01	86	62.0186.01	96	77.0069.00	11
041.637.04	63	044.159.00	67	61.1232.02	86	62.0348.00	12	77.0069.01	11
041.637.05	63	044.170.00	67	61.1267.00	90	62.0655.00	18	77.0095.00	73, 81
041.637.06	63	044.183.00	93	61.1287.00	90	62.0658.00	18	77.0095.01	73, 81
041.642.00	92	044.184.00	93	61.1337.00	91	62.0676.00	85	77.0098.00	73
041.643.00	92	044.184.01	78, 93	61.1373.01	92	62.0685.00	85	77.0098.01	73
041.644.00	92	45.0403.00	60	61.1421.00	85	62.0847.00	85	77.0104.01	59
041.659.01	62	45.0438.00	24	61.1502.00	95	62.0848.00	85	77.0104.03	59
041.660.00	73	47.1147.03	21	61.1503.00	95	65.0208.00	100	77.0121.01	70
041.684.00	62	47.1259.00	20	61.1537.01	95	65.0209.00	100	77.0121.03	70
041.709.00	80	47.2031.00	64	61.1548.00	86	65.0227.00	100	77.0147.00	11
041.806.00	66	048.107.00	, 18	61.1551.00	86	65.0286.02	99	77.0147.01	11
41.1143.00	20	049.001.01	69	61.1553.00	86	65.0286.05	99	77.0254.01	70, 71
41.1144.00	20	049.002.01	69	61.1569.00	86	65.0289.00	100	77.0254.03	70, 71
41.1164.00	20	049.004.01	69	61.1593.00	85	65.0289.01	100	77.0263.00	6
41.1199.00	20	049.005.01	69	61.1632.00	85	65.1625.03	99	77.0263.01	6
41.1299.00	20	049.070.00	6	61.1679.00	85	65.1625.04	99	77.0291.01	59
41.1385.00	20	50.0271.00	3, 24	61.2042.00	96	75.0035.01	73, 74, 75, 76, 77	77.0291.02	59

Index des références

77.0294.00	101, 102	77.1301.00	97	90.1073.00	64	90.1213.00	73	98.0882.00	29
77.0294.01	102	82.0305.04	40	90.1075.00	64, 92	90.1215.00	13	98.0885.00	29
77.0295.00	101	83.0217.00	55	90.1076.00	63	90.1220.00	80	98.0885.12	29
77.0295.01	102	83.0486.00	103, 104	90.1078.00	66	90.1221.00	86	98.1048.00	30
77.0297.00	66	90.0010.00	50	90.1078.01	66	90.1225.00	63, 91	98.1063.00	28
77.0362.01	70	90.0030.00	21	90.1079.00	63, 73	90.1228.00	90	98.1064.00	28
77.0362.03	70	90.0304.00	56	90.1080.00	63	90.1268.00	14	98.1066.00	28
77.0373.00	11	90.0312.00	50	90.1082.00	62	90.1280.00	71	98.1068.00	28
77.0373.01	11	90.0313.00	50	90.1083.00	90	90.1281.00	59	98.1070.00	28
77.0374.00	11	90.0321.00	56	90.1083.01	90	90.1284.00	70	98.1086.00	28
77.0374.01	11	90.0351.00	25	90.1085.00	90	90.1284.01	70	98.1087.00	28
77.0420.01	5	90.0351.01	25	90.1085.01	90	90.1322.00	79	98.1089.00	28
77.0424.00	11	90.0351.02	25	90.1088.00	18	90.1322.01	79	98.1093.00	29
77.0424.01	11	90.0351.03	25	90.1089.00	53	90.1322.02	79	98.1094.00	29
77.0453.00	70, 71	90.0367.01	81	90.1095.00	93	90.1322.03	79	98.1096.00	29
77.0453.01	70, 71	90.0367.02	81	90.1095.01	93	90.1322.04	79	98.1106.00	28
77.0462.00	101	90.0370.00	17	90.1097.00	17	90.1341.00	79	98.1108.00	28
77.0541.00	103	90.0372.00	82	90.1098.00	17	90.1342.00	79	98.1109.00	28
77.0542.00	103	90.0381.00	2, 65	90.1101.00	5	90.1342.01	79	98.1112.00	29
77.0543.00	101	90.0446.00	55	90.1104.00	49	90.1343.00	79	98.1113.00	29
77.0707.00	65	90.0456.00	39, 74, 75	90.1104.01	49	90.1345.00	79	98.1120.00	25
77.0708.00	65	90.0463.01	81	90.1110.00	60	90.1348.00	79	99.0537.01	60
77.0724.00	12, 71	90.0563.00	48	90.1110.01	60	90.1349.00	79	99.0541.01	60
77.0730.01	73	90.0593.00	50	90.1123.00	18	90.1357.00	70, 71	99.0583.00	60
77.0748.00	71	90.1012.00	64	90.1127.00	66	90.1390.00	79	99.0584.00	60
77.0769.00	59	90.1016.00	95	90.1127.01	66	90.1439.00	66	99.0587.00	60
77.0784.00	12	90.1020.00	58	90.1129.00	65, 92	90.1440.00	66	99.0590.00	61
77.0842.00	13	90.1021.00	14	90.1133.00	66	90.1441.00	66	99.0592.00	61
77.0845.00	12	90.1027.03	21	90.1134.00	66	90.1442.00	66	99.0607.00	60
77.0876.00	80	90.1029.00	92	90.1135.00	66	90.1443.00	66	99.0618.00	60
77.0877.00	12	90.1029.01	92	90.1135.01	66	97.0280.02	21	99.0619.00	60
77.0878.00	13	90.1029.02	92	90.1138.00	97	97.0290.00	21	99.0653.01	59
77.0921.00	97	90.1036.00	92	90.1148.00	70	97.0401.01	21	99.0653.03	59
77.0922.00	97	90.1039.00	81, 84	90.1148.01	70	97.0405.00	21	99.0654.02	61
77.0923.00	73	90.1040.00	82	90.1166.00	17	97.0406.00	21	99.0654.03	61
77.0932.00	102	90.1041.00	82	90.1168.00	63	97.0408.01	21	99.0667.00	41
77.0941.01	71	90.1042.00	82	90.1176.00	62	98.0003.00	30	99.0673.01	59
77.0942.00	102	90.1043.00	82	90.1177.00	63	98.0010.00	25, 47	99.0675.01	59
77.0943.00	102	90.1044.00	81	90.1183.00	93	98.0050.01	25	99.0681.00	59
77.0980.00	65	90.1045.00	84	90.1184.00	66	98.0051.00	24	99.0681.01	59
77.0981.00	65	90.1053.00	82	90.1185.00	66	98.0090.01	25	99.0684.00	59
77.1033.00	73	90.1054.00	84	90.1186.00	66	98.0180.00	30	99.0684.01	59
77.1034.00	73	90.1064.00	7	90.1189.00	37	98.0190.00	30	99.0696.01	61
77.1047.00	71	90.1065.00	6	90.1196.00	66	98.0450.00	30	99.0705.00	59
77.1104.00	70	90.1066.00	21	90.1197.00	66	98.0660.00	3	99.0705.01	59
77.1104.01	70	90.1069.00	93	90.1200.00	92	98.0867.02	54	99.0707.01	59
77.1239.00	65	90.1071.00	91	90.1205.00	59	98.0879.00	29	99.0708.01	59
77.1240.00	65	90.1072.00	63, 92	90.1212.00	13	98.0879.12	29	99.0711.00	59

Index des références

99.0711.01	59
99.0719.00	86
99.0724.00	61
99.0738.01	59
99.0739.01	59
99.1135.01	61
99.1136.00	83

Index des références

2-Button Conversion Kit	46
2" Post Adapter	97
2" Post Adapter to 411 Chair	97
2-Way	36
2-Way Valve	78
3-hole tubing	28, 29
3-Way Barbed Plastic Tee	15
3-Way Cartridge	57
4-hole tubing	28, 29
4-Way Barbed Plastic Tee	15
5/16" Water Tubing Barb	15
6-pin fiber-optic	28, 29
8-Button Footswitch	96
8-Function Footswitch	96
12 O'Clock Vacuum Canister Assembly ..	70
300-watt power supply	64
500 Brake Handle Assembly	59
551 Vacuum Tower Assembly	70
6300	82, 83

A

Accessory Holder	59
Accessory Knob Kit	53
Accessory Plate	11
Accessory Power Cord	95
Actuator	58
Actuator Rod Assembly	78
Adapter Insert	98
A-dec 300	12
A-dec 300 Air and Water QD Kit	17
A-dec 300 Cuspidor Bowl	73
A-dec 300 Monitor	13
A-dec 300 Support Center, plug kit	13
A-dec 300 wet/dry foot control	49
A-dec 361 water bottle arm	13
A-dec 500 Air QD Kit	17
A-dec 500 Water QD Kit	17
A-dec 500 Wet/dry foot control	49
A-dec syringe tool	47
Adjustable Screw	57
Air-Actuated	40
Air-Actuated Water Shutoff	40
Air Bleed Filter	57
Air Bleed Restrictor	57
Air Bleed Valve	41
Air Bleed Valve Assembly	78
Air Coolant Flow Adjustment Stem	54

Air Electric Switch	67
Arm Bearing Kit	18
Articulating Headrest	86
Aspirator Tips	69
Assistant's Arm Covers	11
Assistant's Stools	100
Autoclavable Traditional Syringe	44

B

Back Cover	100
Barb	14, 15, 16
Barb/Fittings Service Parts Kit	14
Base Down Stop	94
Blank Whip Cover	11
Bowl Rinse/Cup Fill/Safety Stop	73
Bowl Rinse Spout	74, 76, 77
Bowl Rinse Tube Assembly	77
Bowl Screen	77
Bracket Assembly	81
Bracket Kit	86
Brake Actuator Valve	41
Bulb	1, 2, 3
Bulb Cap	80, 84
Bumper and Plug	103
Button Head Socket Screw	91
Button Valves	46

C

Cable Tie Mount	24
Cable Ties	24
Cap for 2" post on A-dec 300 Monitor ..	13
Cart Height Adjustment Knob	99
Cartridge Needle Valve Kit	57
Cascade	74, 75, 77, 78, 80
Cascade 1040 Chair Adapter	98
Caster	99
Century I control block	37
Century II	56, 57
Century II Control Block	56
Century Plus	54, 55
Century Plus Cartridges	54
Century Plus Control Block	54, 55
Chair Adapter	98
Chair Brake Cable	95
Chair Mount	98
Chair upper structure Hole plug	12
Chip Blower	49, 51

Chip Blower Button	51
Circuit Board	63
Circuit Board Kit	73
Circuit Breaker Kit	82
Clips	24
CM-2.2 module	1
Connector	14
Connector Kit	14
Continental Syringe	45
Contour Stool Back Cover Kit	100
Control Adjustment Stems	53
Control Block Diaphragm	53, 58
Control Block QD Plug	53
Control Valve Assemblies	38
Convolute Tubing	35
Coolant Flow Adjustment Stems	58
Cover	11, 12
Cover, Brake Handle Hole	11
Cup Filler Tube Assembly	77
Cup Holder Assembly	77
Cuspidor	73, 74, 75, 76, 77, 78
Cuspidor Bowl Screens	73
Cylinders	90, 100

D

Data Board	63
Data Cable Adapter	62
Data Cables	63
Decade Chair Mount	98
Deluxe Touchpad	62, 66
Dental Light Intensity Switch	62
Doctor's Holder	59
Doctor's Stools	100
Drive Air Flow Adjustment Stem	54
Dr.'s holder	59
Dry Block Conversion Kit	55
Dry Cartridge	53

E

EA-50	1
EA-51LT	1
EA-52LED	1
EA-53 Replacement Parts	1
Electrical Wiring/Cables	95
End Cap for A-dec 300 Flexarm	12
End Cap For a-dec 361 water bottle arm	13
Exhausting Toggle Valve	37

F

Filter	21, 22
Flexarm End Cap	80
Flex Holder Assembly	59
Floor Box	19
Flow Adjustment Key	53
Flow Adjustment Screw	91
Flow Control Screw	58
Flush Toggle Valve	36
Flush Valve Stem	54
Foot Control	49, 50, 52
Foot Control Cover	52
Foot Control Housing	51
Foot Control I	50, 52
Foot Control III	49, 51
Foot Control III Valve	50
Foot Control III Valve Service Kit	50
Foot Control II Valve	50
Foot Control II Valve Service Kit	50
Foot Control I Valve	50
Foot Control I Valve Service Kit	50
Foot Controls	51
Foot Control Tubing	35, 52

G

Gauge	23
Gauge, Panel Mount Kit	23
Glide Bar	86

H

Handpiece Control Adjustment Keys	54
Handpiece Conversion Kits	57
Handpiece Exhaust Muffler	3
Handpiece Oil Collection System	60
Handpiece Pressure Gauge	3, 24
Headrest Locking Knob	86
Height Adjustment Knob	99
Hex Driver	24
Hex Valve	78
Hex Wrench Set	24
High Volume Evacuator (HVE)	35
Holder	59, 60
Hole Plug	97
Hole Plug for Unit	6
Housing Replacement Kit	81
HVE Holder	59
Hydraulic Oil	91

Index alphabétique

hydraulic Solenoid Kits.....	91	Master Toggle		replacement bulbs.....	2	Touchpads.....	66
I		to Power Supply Kit.....	92	Retaining Clips.....	20	Touchpad Upgrade Kits.....	66
IC Control Block.....	56	Master Toggle Valves.....	37	Retaining Rings.....	52	Touch Up Paint.....	25
ICX® Waterline 2-Liter		Micro Valve.....	51, 78	Rubber Bumper Pad.....	52	Track Light Assembly Tool.....	24
Treatment Tablet.....	6	Monitor Post Cap.....	11	S		Transformer.....	82, 83
ICX® Waterline Treatment Tablet.....	7	O		Saliva Ejector Tubing.....	35	Tray.....	103
Indicator Light Cover.....	12	On/Off Indicator.....	57	Screen.....	72	Tray Holders.....	101, 102, 103
In-Line Check Valve.....	41	On/Off Switch Kit.....	81	Self-Contained Water Cap Assembly.....	7	Tray Pad.....	103
Intensity Switch.....	82, 83, 84	Oral Cups.....	69	Shuttle.....	41	Tri-Flo.....	47
Intermediate Light Post.....	80	O-ring.....	72	Sleeve Clamps.....	16	Tri-Flo Syringe.....	47
J		O-Ring Specification Chart.....	9, 10	Sliding Wedge.....	86	Tri-Flo Syringe Service Kit.....	47
Jumpers.....	62	P		Solenoid.....	62, 63, 73, 91	Trim Ring, telescoping assistant arm...	12
K		Pac 1 Toolkit.....	24	Solenoid Assembly.....	91	Truss Head Slot Screw.....	91
Key Adapter with Screw.....	57	Panel Mount Gauge.....	24	Solids Collector.....	72	Tubing Adapter.....	68
L		PCB.....	62	Spout Gasket.....	73	Tubing Plug.....	16
Lamp Socket Kit.....	81	PCB Standoffs.....	62	Spring Washer Kit.....	55	U	
LED Front Housing.....	79	Performer.....	58	Standard Touchpad.....	66	Umbilical Connectors.....	20
LED Indicator Circuit Board.....	79	Performer Monitor Post Assembly.....	103, 104	Standoffs.....	62	Umbilical hole cover.....	12
LED Light Upgrade.....	79	Phenolic.....	75	Stem Assembly.....	56	Umbilical Tubing.....	35
LED Mode Button.....	79	Pinnacle Screen Kit.....	72	Stools.....	99, 100	Uni-Clamp.....	16
LED Neck Pivot Assembly.....	79	Pins.....	55	Stop.....	95	Uni-Flo Syringe Service Kit.....	47
LED Pivot Shaft Assembly.....	79	Pivot Armrest.....	97	Stop Plate Switch Replacement Kit.....	93	Unitized Holders.....	60
Lever Foot Control Service Kit.....	50	Pivot Stops.....	81	Sub-Miniature Toggle Valve.....	37	Upgrade Kits.....	66
Lift Cylinders.....	90	Plug, Control head touchpad hole.....	12	Surgical Tip Brushes.....	69	Upholstery Kits.....	86
Light Bulb.....	80, 81	Plug kit for a-dec 300 Support Center...	13	Switch Assembly.....	73	Upholstery Screws.....	86
Light Cable Kit.....	82, 84	Plug, Monitor cable hole in light post...	12	Switch Service Kit.....	82	V	
Light Circuit Board.....	83	Ported 3-Way.....	40	Swivel Brake Kit.....	95	Vacuum Canister Assembly.....	71
Light Flexarm end cap.....	12	Positioning Sensor Kit.....	93	Syringe Flow Control Barb.....	54	Vacuum canister blank cover.....	12
Light Handle Kit.....	79, 81	Post Trim Rings.....	73	Syringe Kits.....	44	Vacuum Tubing.....	35
Light post cover Plugs.....	80	Potentiometer.....	94	Syringe Service Kit.....	46	Vacuum Valve Brushes.....	69
Light Post Trim Ring.....	81	Power Supply.....	65, 67	Syringe Tool.....	25	Valve Core Service Kit.....	47
Light Shield.....	80, 81, 84	Power Supply cables.....	65	Syringe Valve Kit.....	46	Valve Cover Kit.....	56
Light Shield And Gasket Kit.....	79	Pressure Indicator.....	42	T		Valve Kit.....	37
Light Yoke Plug.....	82	Priority Chair Mount.....	98	Tee Check.....	41	Valve Test Syringe.....	25
Limit Switch.....	78, 93	P-trap.....	73	Tension Springs.....	55	VE-10 tubing.....	1
Liners.....	69	Push-Button Valve.....	36	Terminal Barb Assembly.....	44, 45	Vertical adjustment Kit.....	80
Loctite®.....	25	Push On Knob Kits.....	56	Tie Bolt Kit.....	55, 56	Vitreous China.....	74
Lubricant.....	25	Q		Tilt Actuators.....	90	W	
Lubriplate.....	83	Quad Voltage Intraoral Light Source.....	2	Tilt Cylinders.....	90	Wago Connector.....	62
M		Quick disconnect Kit.....	97	Toggle.....	42	Washers.....	84
Manifold.....	62	R		Toggle Switch Kit.....	84	Water Bottle Cap.....	6
Manual Shutoff.....	40	Radius Flat-Panel Mount.....	103	Toggle Valve.....	36, 37, 78	Water Bottle Tubing Grommet.....	6
		Relay Module.....	63	Tooth Dryer.....	2	Water Cartridge.....	53

Index alphabétique

Water Coolant Valve	37
Water Retraction Tester	24
Water Valve Cover	56
Wet/dry foot control iii	49
Wet/Dry Toggle	51
Whip Assembly Kit	59
Whip Components.	59
Whip hook	59



A-dec Headquarters

2601 Crestview Drive
Newberg, Oregon 97132
United States

Tel: 1.800.547.1883 within USA/CAN
Tel: +1.503.538.7478 outside USA/CAN
Fax: 1.503.538.0276

a-dec.com
a-dec.biz

A-dec Australia

Unit 8, 5-9 Ricketty Street
Mascot, NSW 2020
Australia

Tel: 1.800.225.010 within AUS
Tel: +61 (0)2 8332 4000 outside AUS
Fax: +61.(0)2.9699.4700

A-dec China

A-dec (Hangzhou) Dental Equipment Co., Ltd.
528 Shunfeng Road, Qianjiang Economic Development Zone
Hangzhou 311106
Zhejiang, China

Tel: +1.503.538.7478

A-dec United Kingdom

Austin House, 11 Liberty Way
Nuneaton, Warwickshire CV11 6RZ
England

Tel: 0800.ADEC.UK (2332.85) within UK
Tel: +44.(0)24.7635.0901 outside UK
Fax: +44.(0)24.7634.5106
